KRUUSE FIBREGLASS CASTING TAPE

White and Blue

GB KRUUSE Fibreglass Casting Tape

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
Instructions for use:

Usage and composition:

KRUUSE orthopaedic casting tape is designed for use in the construction of the most common orthopaedic casts, as well as prosthetics and orthotic devices. Specific and suitable applications should be the responsibility of onsite qualified medical professionals. KRUUSE orthopaedic casting tape consists of a specially knitted fibreglass fabric impregnated with specially developed water-activated polyurethane resin. Exposure to moisture or water initiates a chemical reaction that causes the material to be hardened. Finished casts produced from KRUUSE orthopaedic casting tape are lightweight, durable, strong, porous and radiolucent.

Precautions:

1. Protective or surgical gloves must be worn while handling the tape, as the polyurethane resin in KRUUSE orthopaedic casting tape will firmly adhere to unprotected skin and clothing. Care should be taken to avoid coming into contact with unprotected areas of the patient's skin during application. Applied tape should be allowed to fully cure for about 30 minutes before touching unprotected skin.
2. In the event resin comes into contact with skin or instruments, remove it immediately using acetone or alcohol before the resin sets.
3. As KRUUSE orthopaedic casting tape is slimline and lightweight, patients should be advised to avoid activities that may affect the integrity of the cast and the healing process.
4. A certain amount of joint stiffness is to be expected in certain cases, such as long-term immobilisation in a cured cast.
5. Increased heat causing discomfort may result from excessive layers of tape, dipping into water above 80 °F (27 °C), squeezing too much water from the tape before application, using an elastic bandage on a freshly applied cast for more than one minute, covering the cast with a blanket or otherwise restricting airflow around the cast. (Excessive heat may cause burns.)
6. Respiratory sensitivity may occur. It is recommended that casting tape is applied in a ventilated area.

Application of KRUUSE orthopaedic casting tape:

1. Apply one or two layers of stockinette and cast padding over the affected part. If desired, additional cast padding may be used for greater protection of pressure points and bony prominences. In case of immersion in water, use of non-absorbent undercast material or a water-resistant cast cover after the application is recommended.
2. Protective gloves must be worn while wrapping the casting tape to avoid the resin in the tape transferring and adhering to the skin.
3. Select the desired width of the KRUUSE orthopaedic casting tape. Open only one roll at a time, room humidity may initiate hardening.
4. Immerse the roll of tape in room temperature water 68-77°F (20-25°C). **Note:** Squeezing the roll 3-4 times while it is submerged gives a set time of approximately 3 minutes. To lengthen the set time, immerse the roll and remove it immediately without squeezing. This technique gives a set time of approximately 4 minutes. If the roll is not squeezed while submerged in water, warmer water will shorten the set time while cooler water will lengthen it. If the roll is immersed in water above 80°F (27°C), an increase in heat may be noted. A non-immersed dry cast is recommended for a longer set time.
5. Construct the cast by wrapping in spiral motions, overlapping the previous layer by one-half to two-thirds the width of the tape. As with any cast, take care to avoid excessive tightness. To ensure a snug fit, wrap with slight tension on the tape.
6. Just three or four layers provide a strong, non-weight-bearing cast. Additional layers are generally not necessary for strength and will reduce porosity.
7. Smooth and rub the tape for good lamination between layers. Soon after smoothing, the cast can be moulded.
8. About 3-4 minutes after immersion in water, KRUUSE orthopaedic casting tape is sufficiently rigid to prevent further moulding.
9. About 20 minutes after setting, the cast is able to bear weight.
10. Walking heels, cast shoes and bracing devices can all be used with KRUUSE orthopaedic casting tape. The cast can be modified by reinforcing or patching with

additional tape for therapeutic reasons. The fresh tape should be overlapped on itself to secure in place.

Shelf life and storage requirements:

1. Each roll of tape is packaged in a vacuum-sealed or dried air-injected foil pouch. Care should be taken to avoid puncturing the pouch as this will cause the tape to harden prematurely in the package. If the packaging is punctured due to incorrect storage or handling and it feels hard, dispose of it and do not use the tape. It is recommended squeezing the packaging gently prior to opening to check it is suitable for application.
2. If the packaging becomes wet, it must be dried carefully with a towel and then stored in a dry cool place. For best results, store the product between 59-86 °F (15-30 °C). For shelf life, refer to the expiration date on the packaging.
3. Each box and pouch is marked with a lot number and expiration date. For optimum performance, rotating stock is recommended.

Cast removal and disposal:

To remove the cast, use an oscillation-type cast saw with a vacuum. While the dust generated is less than that generated by plaster of Paris casts, take precautions to avoid breathing in the dust, as this may cause respiratory irritation or sensitivity. The KRUUSE orthopaedic cast can have openings inserted into it, be trimmed or also removed using traditional saw equipment. Ensure that the blade is always sharp to avoid any heat build-up. Cast edges are abrasive, so take care to avoid scratches. Casts may be incinerated or disposed of to landfill in accordance with local regulations.

Flammability:

KRUUSE orthopaedic casting tape meets class I flammability requirements of title 16 CFR 1610 (CS191-53).

Instructions for cast patients:

1. Patient should be advised to avoid activities that may be too vigorous or which may interfere with the reduction of the fracture or the integrity of the cast.
2. If the cast becomes soiled, it can be cleaned with a damp cloth. If there is a rough edge on the cast, it can be smoothed off with a nail file or emery board.

3. Immersion in water is not recommended. In case of immersion, instruct the patient on proper drying techniques. To dry a wet cast, blot with a towel to absorb excess moisture, then use a hair dryer until the padding and undercast material is fully dry.
4. Care should be taken to prevent the accumulation of foreign materials, such as sand or dirt under the cast, as they may cause skin irritation.

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