KRUUSE

KRUUSE VET-10 URINE STRIPS

Urobilinogen

Glucose

Bilirubin

Ketones

S.G

Blood

pH

Protein

Nitrite

Leukocytes

100/pk

290798

KRUUSE VET-10 URINE STRIPS

Reagent Strips for Veterinary Urinalysis

INTENDED USE

Urinalysis is an important tool in disease detection, as well as monitoring and screening animal health. Abnormalities can be indicative of diseases of the urinary system as well as other organ systems. Complete urinalysis involves both macroscopic and microscopic assessment. This is typically performed by gross visual assessment of the urine, microscopic examination and chemical evaluation. Several chemical parameters can be measured and the KRUUSE VET-10 Urine Strips include the following parameters: bilirubin, blood, glucose, ketones, pH, protein, specific gravity, urobilinogen, leukocytes and nitrite. For veterinary use only.

STORAGE AND HANDLING

Store in a cool, dry place at temperatures between 2°C~30°C. Do not store the strips in a refrigerator. Store away from moisture and light. When stored in the original container, the product is stable up to the expiry date printed on the box and container. Replace the cap immediately and tightly after removing test strips, and do not remove desiccant from container. Do not touch test areas of urine reagent strips. Do not open container until ready to perform the test. Discoloration or darkening of the test pads may indicate deterioration. If this is case, or if test results are questionable or inconsistent with expected finding, confirm that the product is within its expiration date and additionally check if the test strip is reacting properly, using known negative and positive control materials. Do not use the strips after the expiry date. Note once the container has been opened, the remaining strips are stable for up to 6 months.

QUALITY CONTROL

For best results, performance of reagent strips should be confirmed by testing known negative and positive specimen or controls when a new container is first opened. Each laboratory should establish its own standards of performance.

SAMPLE COLLECTION

Urine should be collected in a clean, dry container, free of any disinfecting or cleaning chemicals. Collect the urine sample using the protocol of the clinic or hospital. Be aware that samples collected by normal urination may be contaminated from surrounding environment and that samples collected by cystocentesis may contain red blood cells.

SAMPLE HANDLING

In order to obtain accurate results, the urine collection, storage and handling must be sterile and follow standard procedures. The dipstick analysis should be performed as soon after collection as possible (ideally within 30 minutes of collection) and the sample should be well mixed prior to testing. If, for some reason, the test cannot be performed immediately, the sample may be covered and refrigerated. The sample should return to room temperature prior to testing.

TESTING METHODS

Remove the dipstick from the container. It is important not to touch the reagent areas of the strip as this may alter test results. Each reagent area should be immersed in urine by dipping. Remove excessive urine to prevent dilution of reagents or mixing of reagents between pads. This can be achieved by tilting the strip and allowing the urine to run off the edges. While blotting excessive urine, ensure the chemicals from the different parameters do not mix. The reagent pads should be read at the specified times. These times are different for each test. Compare the blocks to the corresponding colour chart provided. Urine discoloration may create difficulty in visually interpreting the test results. If the urine is noticeably discolored, the sample may be centrifuged and the supernatant used for analysis.

UROBILINOGEN

Urobilinogen is formed by intestinal bacteria from the breakdown of conjugated bilirubin. Urobilinogen is usually excreted in feces, however a small amount may be reabsorbed and excreted in urine. This test is not of significant value in animals. High concentrations of biliary pigments may occur in hemolytic crisis, or cases of hepatic or intestinal dysfunction. A false negative test result may occur if the sample is old, as urobilinogen is very unstable when exposed to light and air.

GLUCOSE

Glucose is not detectable in the urine of healthy dogs or cats, as glucose passes freely through the glomerular filter and is resorbed by the proximal tubules. If glucosuria is present, it is due to either an excess amount of glucose reaching the tubules that cannot be resorbed or, less commonly, decreased tubular resorptive function. Glucosuria may be either persistent or transient and multiple tests may be needed for differentiation of these conditions. Persistent causes of glucosuria include: diabetes mellitus, administration of glucose containing fluids, hyperadrenocorticism, hyperpituitarism or acromegaly. Other diseases that may result in transient hyperglycemias leading to glucosuria include: hyperthyroidism, acute pancreatitis, stress (especially in cats), postprandial, and administration of certain drugs. Moderately high concentrations of ketones may cause false negative test results if the amount of glucose is only slightly elevated. Cold urine (refrigerated specimens) or expired reagent strips may also result in false negative test results.

BILIRUBIN

Bilirubin is produced from the breakdown of hemoglobin, transported to the liver bound to albumin, and conjugated with carbohydrates by hepatocytes. Only conjugated bilirubin is found in urine. Excess bilirubin may be produced when red blood cells are destroyed, or in liver disease, including bile duct obstruction. Conjugated bilirubin is detected in urine if the renal threshold is exceeded. The renal threshold in dogs, especially males, is lower than that of other species. Bilirubin is very unstable when exposed to room air and light. Thus, urine specimens should be tested soon after collection. Positive test results may be observed in concentrated urine of healthy dogs. In dogs, the renal threshold for bilirubin is low and renal tubules are able to break down heme and produce some renal bilirubin, therefore slight bilirubinuria can be a normal finding in dogs with concentrated urine. However, bilirubinuria is always abnormal in cats. Bilirubinuria may indicate: liver disease, bile duct obstruction, starvation, hemolysis, or pyrexia.

KETONES

Acetone, acetoacetic acid, and beta-hydroxybutyric acid are ketones. Glomeruli freely filter ketones and the tubules then resorb them completely. If the tubular resorptive capacity is saturated, then the ketones are incompletely resorbed, resulting in ketonuria. Ketonuria occurs quickly in younger animals and is more easily detected than ketonemia. Ketonuria does not signify renal disease, but rather excessive lipid or defective carbohydrate metabolism. Dipstick tests are semiquantitative and only detect acetone and acetoacetic acid. Ketonuria may be caused by starvation, insulinoma, diabetic ketoacidosis, persistent hypoglycemia, high fat low carbohydrate diets, and glycogen storage disease.

SPECIFIC GRAVITY

Urine specific gravity is based on the ratio of weight of urine to weight of an equivalent volume of pure water. This test is used to measure tubular function. Although dipstick strips do have a method of approximating specific gravity, this measurement is best made with a refractometer.

BLOOD

The test pad will react positively in the presence of red blood cells, free hemoglobin, or free myoglobin. Hemoglobin usually is bound and is too large to pass through the glomerular filter. If the renal threshold is exceeded, the hemoglobin can pass into the urine. Myoglobin is not bound and passes freely through the glomerular filter. Myoglobin can be detected in urine before a change in plasma color is apparent. The presence of free red blood cells results in a positive test when blood cells lyse, and hemoglobin is released. Healthy animals should have negative test results. A positive test indicates hematuria, hemoglobinuria or myoglobinuria. Most commonly, hematuria is the cause of the positive test result while myoglobinuria is rare. Hematuria can be caused by trauma, infection, inflammation, infarction, calculi, neoplasia, or a coagulopathy anywhere along the urinary

tract. In cases of hematuria, the urine is red and cloudy, but will clear if centrifuged. Microscopic evaluation of the urine sediment will reveal red blood cells. Hemoglobinuria, on the other hand, will have reddish brown urine that does not become clear after centrifugation. The microscopic evaluation of urine sediment will not reveal red blood cells.

PH

The pH of urine can vary depending on an animal's diet as well as its acid-base status. Animals that primarily eat high protein meat-based diets will have acidic urine, as animals that eat more vegetable-based diets will have alkaline urine. The urine sample should be fresh as the urine becomes more alkaline on standing due to the conversion of urea to ammonia by bacteria (if present), and loss of CO₂. Causes of acidic urine include: a meat diet, systemic acidosis, hypochloremia, and administration of acidifying agents. Urine with high concentrations of glucose may have a lower pH, due to the bacterial metabolism of glucose and the production of ammonia which lowers

pH. Causes of alkaline urine include: vegetable-based diet, a bacterial infection of urease-producing bacteria, systemic alkalosis, urine exposed to room air for an extended time (loss of CO₂), and administration of alkalinizing agents including citrate or NaHCO₃. Urine pH also may provide a good predictive assessment of crystal and stone morphology as certain crystals and stones form in either acidic or alkaline environments. Uric acid, cystine, and calcium oxalate crystals are found in acidic urine, as struvite, calcium carbonate, calcium phosphate, ammonium biurate, and amorphous phosphate crystals are found in alkaline urine.

PROTEIN

Dogs and cats normally have small proteins that pass through the glomerular filter, however, a majority of these proteins are resorbed by the renal tubules. Only a very small amount of protein is normally excreted in the urine, which is not usually clinically detectable. Proteinuria may be due to hemorrhage, infection, intravascular hemolysis, or renal disease. Hemorrhage is confirmed by a positive occult blood reaction on the dipstick and the presence of red blood cells in the sediment. A urinary infection or cystitis can be confirmed by observing bacteria and white blood cells on sediment examination. Cases of intravascular hemolysis have hemoglobinuria leading to a positive occult blood test. Proteinuria of renal disease may be due to glomerular and/or tubular lesions. If the proteinuria is due to renal disease, the occult blood test will be negative, and the sediment may or may not contain casts. Determination of the urine protein/urine creatinine ratio is helpful in confirming renal proteinuria. Protein results must be analyzed with the urine-specific gravity. Trace proteinuria may represent significant protein loss with low specific gravity, but not with high specific gravity. False positive protein reactions may occur with alkaline urine. Samples containing urease-producing bacteria may have an elevated pH resulting in a false positive test result.

NITRITE

The nitrite portion of the dipstick analysis has limited value in veterinary medicine. Nitrites occur in urine during some bacterial infections. In order to achieve an accurate positive test result, the urine must have been retained in the bladder for at least 4 hours. Therefore, it is best to collect a (first) morning sample or ensure the patient has not urinated in at least 4 hours. A positive test indicates a bacterial infection. Gram-negative rods are more likely to produce a positive test response. Negative test results do not exclude infection. The urinary tract infection may involve organisms that do not convert nitrites or the urine may not have been held in the bladder greater than 4 hours.

LEUKOCYTES

The leukocyte test detects the presence of white blood cells or partial cells in the urine. In dogs, this test is indicative of pyuria but false negative test results often occur. False positive test results often occur in cats, and this test is clinically unreliable. False positive test results also may occur in the event of fecal contamination. False negative test results may develop if the patient has been treated with high doses of tetracycline or other antibiotics. Glucosuria or increased urine-specific gravity may cause false negative test results. False negative test results may be observed with voided urine samples obtained from animals with pyometra or prostatitis.

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