

KRUUSE Postoperative Care



KRUUSE Postoperative Care

This chapter contains our range of postoperative care products used by veterinarians and technicians worldwide.

We offer a wide range of disposable collars for cats and dogs in the original BUSTER design and other products to prevent the animal from reaching certain areas of the body.

BUSTER Collars

- BUSTER Quick Collar
- BUSTER Complete Collar
- BUSTER Classic Collar
- BUSTER Classic Transparent Collar
- BUSTER Comfort Collar
- BUSTER Clic Collar
- BUSTER Design Collar
- BUSTER Dog Collar for Brachycephalic Breeds
- Inflatable Collars



BUSTER Quick Collar

The BUSTER Quick Collar has a touch fastener system for easy customised adaptation to each animal. The collar is padded with a soft outer silicone edge to avoid marking or scraping and a soft inner silicone edge for a comfortable fit. The collar is provided with neck holes for secure fastening to the dog's neck collar with ie. gauze or a strip.

- Quick custom fit for the individual animal
- Easy application and removal
- Full vision
- Full body coverage
- Made of soft silicone and 100% polypropylene
- 7 sizes and a Starter Pack available

Recommended use for BUSTER Quick Collar

- Postoperative care
- Wound care
- Spay/neuter recovery
- Hot spots and other skin ailments
- Skin diseases



Unique features and benefits

Safety and perfect fit

The collar is provided with neck holes for secure fastening to the dog's neck collar with ie. gauze or a strip.

Easy application

The collar has a touch fastener system for easy and customized adaptation to each animal.

Soft padded edges

The collar is padded with a soft outer silicone edge to avoid scratching humans and furniture and a soft inner silicone edge for optimal a comfortable fit.


Full body coverage

The collar offers a full body coverage and protection preventing the animal from reaching certain areas of the body.

Instruction for use


Step 1

- Fold the collar around the neck
- Mark the selected position
- Leave two fingers width between the collar and neck
- Keep the touch fastener mated throughout the sizing procedure



Step 2

- Remove the adhesive liner from the touch fastener




Step 3

- Apply the adhesive part to the marked position
- Cut off redundant collar material



Step 4

- Place the collar on the pet's neck



Instruction for use
Scan QR code to
watch the video



BUSTER Quick Collar

Caution to Pet Owners

It is normal that pets are responding negatively the first time, they are introduced to the e-collar. To ensure your pet's safety while wearing the e-collar keep it indoors in a restricted area. If you are concerned about your pet's behavior or recovery contact your veterinarian.

Collar Maintenance

Wipe with a damp cloth.



Cat. No	Size cm / in	Pk. Size	Suitable for	Kg	lb
273511	7.5 cm / 3"	5/pk	Pomeranian, Chihuahua	1-3	2-6
273512	10 cm / 4"	5/pk	Yorkshire terrier, Bichon Frise, Maltese	3-5	6-11
273513	12.5 cm / 5"	5/pk	Danish/Swedish Farmdog, Papillon, Cairn Terrier, Border Terrier	5-10	11-22
273514	15 cm / 6"	5/pk	Dutch Kooikerhondje, Scottish Terrier, Beagle, Irish Terrier	10-15	22-33
273515	20 cm / 8"	5/pk	Australian Shepherd, English Springer Spaniel, Staffordshire Bull Terrier	15-25	33-55
273517	25 cm / 10"	5/pk	Labrador Retriever, Siberian Husky, Irish Setter	25-30	55-66
273518	30 cm / 12"	5/pk	Golden Retriever, Boxer, Collie, German Shepherd	25-35	55-77

273519 Starter Pack of size 7.5 - 30 cm 7/pk

BUSTER Complete Collar

Superior comfort and protection for post-op and more.

The BUSTER Complete Collar is an ideal choice for patients recovering from surgery or minor injuries. It is also suitable to protect lesions resulting from or aggravated by excessive licking, biting or chewing.

Features and benefits

- Quick and easy to custom fit
- Easy application and removal
- Full vision
- Full body coverage
- Made of soft perforated neoprene and 100% polypropylene
- Available in 7 sizes and a starter pack

Recommended uses

- Postoperative care
- Wound care
- Spay/neuter recovery
- Hot spots and other skin ailments
- Skin diseases



Scan QR code to
watch the video



BUSTER



Cat. No	Size	Pk Size	Suitable for	Kg
273551	7.5 cm / 3"	5/pk	Pomeranian, Chihuahua	1-3 kg
273552	10 cm / 4"	5/pk	Yorkshire Terrier, Bichon Frise, Maltese	3-5 kg
273553	12.5 cm / 5"	5/pk	Danish/Swedish Farmdog, Cairn Terrier, Border Terrier	5-10 kg
273554	15 cm / 6"	5/pk	Dutch Kooikerhondje, Scottish Terrier, Beagle, Irish Terrier	10-15 kg
273555	20 cm / 8"	5/pk	Australian Shepherd, English Springer Spaniel	16-25 kg
273556	25 cm / 10"	5/pk	Labrador Retriever, Siberian Husky, Irish Setter	25-30 kg
273557	30 cm / 12"	5/pk	Golden Retriever, Boxer, Collie, German Shepherd	25-35 kg
273558		7/pk	Starter Pack of sizes 7.5 - 30 cm	

BUSTER Clic Collar

The quick and simple clic closing system makes the application of the collar as easy as buttoning your shirt. Fully transparent material will allow the animal all-round vision. Soft and flexible.

Made of environmentally friendly material. 100% PP - polypropylene.

Cat. No	Description
273480	BUSTER Clic Collar, 7.5 cm, 10/pk
273481	BUSTER Clic Collar, 10 cm, 10/pk
273482	BUSTER Clic Collar, 12.5 cm, 10/pk
273483	BUSTER Clic Collar, 15 cm, 10/pk
273484	BUSTER Clic Collar, 20 cm, 10/pk
273485	BUSTER Clic Collar, 25 cm, 10/pk
273486	BUSTER Clic Collar, 30 cm, 10/pk
273488	BUSTER Clic Collar, 35 cm
273487	BUSTER Clic Collar, 40 cm
273490	BUSTER Clic Collar combi pack, 20/pk 2 x 7.5 cm, 3 x 10 cm, 3 x 12.5 cm, 3 x 15 cm, 3 x 20 cm, 3 x 25 cm and 3 x 30 cm
273489	BUSTER Clic Collar, assorted sizes from 7.5 to 30 cm, 7/pk



BUSTER Comfort Collar

Collar with rubber edge - equipped with a new quick-fastening system. The collar is transparent which enables the animal full vision. Soft and flexible.

Made of environmentally friendly material. 100% PP - polypropylene.

Cat. No	Description
273900	BUSTER Comfort Collar, 7.5 cm, 5/pk
273901	BUSTER Comfort Collar, 10 cm, 5/pk
273902	BUSTER Comfort Collar, 12.5 cm, 5/pk
273903	BUSTER Comfort Collar, 15 cm, 5/pk
273904	BUSTER Comfort Collar, 20 cm, 5/pk
273905	BUSTER Comfort Collar, 25 cm, 5/pk
273906	BUSTER Comfort Collar, 30 cm, 5/pk
273907	BUSTER Comfort Collar, assorted sizes from 7.5 to 30 cm, 7/pk





Illustration of BUSTER Clic / Comfort Collar



How to assemble the BUSTER Clic / Comfort Collar

BUSTER Classic Collar, Transparent

Our well-known classic dog collar is also available in a transparent version. Flexible and easy to handle. Easy to fit on the patient. Available in 9 sizes to ensure perfect fit for any dog.

Made of environmentally friendly material. 100% PP - polypropylene.

Cat. No	Description
273391	BUSTER Transparent Dog Collar, 7.5 cm, 10/pk
273392	BUSTER Transparent Dog Collar, 10 cm, 10/pk
273393	BUSTER Transparent Dog Collar, 12.5 cm, 10/pk
273394	BUSTER Transparent Dog Collar, 15 cm, 10/pk
273395	BUSTER Transparent Dog Collar, 20 cm, 10/pk
273396	BUSTER Transparent Dog Collar, 25 cm, 10/pk
273397	BUSTER Transparent Dog Collar, 30 cm, 10/pk
273399	BUSTER Transparent Dog Collar, 35 cm
273398	BUSTER Transparent Dog Collar, 40 cm
273389	BUSTER Transparent Dog Collar, assorted sizes from 7.5 to 30 cm, 7/pk



BUSTER Classic Collar

Protective disposable collars for dogs.

An inexpensive solution to prevent the patient from reaching and damaging wounds on head and body. Flexible and easy to handle. Easy to fit on the patient. Available in 9 sizes to ensure perfect fit for any dog.

Made of environmentally friendly material. 100% PP - polypropylene.

Cat. No	Description
273310	BUSTER Classic Dog Collar, 7.5 cm, 10/pk
273320	BUSTER Classic Dog Collar, 10 cm, 10/pk
273325	BUSTER Classic Dog Collar, 12.5 cm, 10/pk
273330	BUSTER Classic Dog Collar, 15 cm, 10/pk
273340	BUSTER Classic Dog Collar, 20 cm, 10/pk
273350	BUSTER Classic Dog Collar, 25 cm, 10/pk
273360	BUSTER Classic Dog Collar, 30 cm, 10/pk
273368	BUSTER Classic Dog Collar, 35 cm
273367	BUSTER Classic Dog Collar, 40 cm
273369	BUSTER Classic Dog Collar, assorted sizes from 7.5 to 30 cm, 7/pk





Illustration of BUSTER Classic Collar




How to assemble the BUSTER Classic Collar

BUSTER Design Collar Bones

Make the sending home of the animal more fun with the new BUSTER Design Collar with printed bones. Watch the surprise and happiness of your clients, when the dog wears printed collar in eye-catching colours.

- Quick BUSTER Clic Fastening system
- Color-coded according to size
- Full body protection
- Multiple adjustment points for proper fit
- Available in 7 sizes



Recommended use for BUSTER Design Collar Bones


- Postoperative care
- Wound care
- Spay/neuter recovery
- Hot spots and other skin ailments
- Skin diseases

Cat. No	Size	Pk. Size	Colour
273731	7.5 cm	10/pk	
273732	10cm	10/pk	
273733	12.5 cm	10/pk	
273734	15 cm	10/pk	
273735	20 cm	10/pk	
273736	25 cm	10/pk	
273737	30 cm	10/pk	

273538 BUSTER Design Collar Bones, assorted sizes from 7.5 to 30 cm, 7/pk



Instructions for use
Scan QR code to watch the video



Size	Length in cm / Inch	Opening in cm / Inch	Neck-hole in cm / Inch	Suitable for	Kg
7.5 cm	7.5 cm / 2.95"	13-15 cm / 5.1"-5.9"	7-8 cm / 2.7"-3.1"	Pomeranian, Chihuahua	1-3 kg
10 cm	10 cm / 3.93"	15-18 cm / 5.9"-7.0"	7-8 cm / 2.7"-3.1"	Yorkshire terrier, Bichon Frise, Maltese	3-5 kg
12.5 cm	12.5 cm / 4.92"	19-21 cm / 7.4"-8.2"	9-10.5 cm / 3.5"-4.1"	Danish/Swedish farmdog, Papillon, Cairn Terrier, Border Terrier	5-10 kg
15 cm	15 cm / 5.90"	20-22 cm / 7.8"-8.6"	10-12 cm / 3.9"-4.7"	Dutch Kooikerhondje, Scottish Terrier, Beagle, Irish Terrier	10-15 kg
20 cm	20 cm / 8.87"	23-27 cm / 9.0"-10.6"	12-14 cm / 4.7"-5.5"	Australian Shepherd, English Springer Spaniel, Staffordshire Bull Terrier	15-25 kg
25 cm	25 cm / 9.84"	34-42 cm / 13.3"-16.5"	14-16 cm / 5.5"-6.2"	Labrador Retriever, Siberian Husky, Irish Setter	25-30 kg
30 cm	30 cm / 11.81"	41-48 cm / 16.1"-18.8"	15-18 cm / 5.9"-7.0"	Golden Retriever, Boxer, Collie, German Shepherd	25-35 kg

BUSTER Design Collar Under Construction

The BUSTER Design Collar is an eye-catching printed collar with a roadworks theme and road signs with cute statements wishing the dog good luck.

The BUSTER Design Collar is transparent and has the quick and simple clic closing system, which is well known from the BUSTER Clic Collar.

Made of environmentally friendly material. 100% PP - polypropylene.

Cat. No	Description
273520	BUSTER Design Dog Collar, 7.5 cm "Under Construction", 10/pk
273521	BUSTER Design Dog Collar, 10 cm "Under Construction", 10/pk
273522	BUSTER Design Dog Collar, 12.5 cm "Under Construction", 10/pk
273523	BUSTER Design Dog Collar, 15 cm "Under Construction", 10/pk
273524	BUSTER Design Dog Collar, 20 cm "Under Construction", 10/pk
273525	BUSTER Design Dog Collar, 25 cm "Under Construction", 10/pk
273526	BUSTER Design Dog Collar, 30 cm "Under Construction", 10/pk
273527	BUSTER Design Dog Collar, "Under Construction", assorted sizes from 7.5 to 30 cm, 7/pk



BUSTER Dog Collar for Brachycephalic Breeds

Our range of BUSTER Collars has been extended with the new BUSTER Dog Collar for Brachycephalic Breeds. Finally, you can get a suitable collar for the brachycephalic dogs without having to adjust and cut the collar.

- Quick fastening system
- Full body coverage
- Multiple adjustment points for proper fit

Recommended use for BUSTER Dog Collar for Brachycephalic Breeds

- Postoperative care
- Wound care
- Spay/neuter recovery
- Hot spots and other skin ailments
- Skin diseases


The BUSTER Dog Collar for Brachycephalic Breeds is transparent and has the quick and simple "Clic" closing system which is well-known from the BUSTER Clic Collar.

The collars are made of 100% polypropylene.



Cat. No	Size	Pk. Size	Suitable for	Kg
273500	XS	10/pk	Pekingese, Japanese Chin, Brussels griffon, Chihuahua	2-6 kg
273501	S	10/pk	Pug, Boston terrier, Cavalier King Charles spaniel, Shih Tzu	6-9 kg
273502	M	10/pk	French Bulldog, Lhasa Apso	9-16 kg
273503	L	10/pk	English Bulldog, American bulldog, Boxer, Bullmastiff	18-25 kg

273504 BUSTER Collar for Brachycephalic Breeds, assorted sizes, 2 x XS, 3 x S, 3 x M, 2 x L



Size	Neck-hole cm	Neck-hole inches	Opening cm	Opening inches	Length cm	Length inches
XS	8-9 cm	3.1-3.5"	17-20 cm	6.6-7.8"	10 cm	3.9"
S	9.5-10.5 cm	3.7-4.1"	23-25 cm	9.0-9.8"	12 cm	4.7"
M	12-14 cm	4.7-5.5"	30-33 cm	11.8-12.9"	16 cm	6.2"
L	14-17 cm	5.5-6.6"	36-40 cm	14.1-15.7"	20 cm	7.8"

Scan QR code to watch the video



Collar Rack

Stainless steel rack for storing and displaying the BUSTER collars in the clinic. Made for seven sizes.

Cat. No	Description
273370	Stainless steel rack for collars



Strap to Fix Collars

To fix BUSTER collars or to be used as an identity band.

Cat. No	Description
273300	Strap to fix BUSTER collars, 10/pk



BUSTER Neckstrap

Practical and flexible neckstrap which is used in connection with a BUSTER collar when the dog's normally worn collar is not suitable for this purpose.

Cat. No	Description
273492	Neckstrap for 7.5 and 10 cm, BUSTER Clic Collar, 10/pk
273494	Neckstrap for 12.5 and 15 cm, BUSTER Clic Collar, 10/pk
273497	Neckstrap 20, 25 and 30 cm, BUSTER Clic Collar, 10/pk
273498	Neckstrap for 40 cm, BUSTER Clic Collar, 10/pk



BUSTER Soft Collar

For cats and small dogs. Comfortable, soft, non-woven material which is waterrepellent and easy to clean. Particularly suitable for cats due to the shape and flexibility of the collar.

Cat. No	Description
273384	BUSTER Soft Collar for cats and small dogs, 10/pk



BUSTER Rabbit Collar

Protective collar for rabbits. Made of flexible and strong polyethylene. An inexpensive solution to prevent the animals from reaching wounds on head and body.

Cat. No	Description
273390	BUSTER Rabbit Collar, 10/pk



BUSTER Bird Collar

Disposable collar for birds. Prevents the bird from picking feathers and wounds. Easy to fit on the bird with an office stapler. Transparent material.

Cat. No	Description
273375	BUSTER Bird Collar, 80 mm, 10/pk
273376	BUSTER Bird Collar, 100 mm, 10/pk
273377	BUSTER Bird Collar, 140 mm, 10/pk
273378	BUSTER Bird Collar, 200 mm, 10/pk
273372	BUSTER Bird Collar, combi, 15/pk
273386	BUSTER Parrot Collar, normal, 35 mm
273387	BUSTER Parrot Collar, L, 45 mm

273386-87



273372-78



Inflatable Collars

Soft and comfortable collar to prevent the dog from reaching injuries, stitches, rashes or wounds. The protective inflatable collar has an inflatable inner pouch that is inflated through an air valve. The inner pouch is protected by an outer nylon or PVC cover to resist bites and scratches.

Features

- Will not inhibit the animal's field of vision
- Durable long-term use
- Comfortable for the animal
- Protective against injuries, rashes and post-surgery procedures
- Will not mark or scrape furniture
- Scratch and bite-resistant



BUSTER Inflatable Collar Nylon

Cat. No	Description	External diameter	Internal diameter	Suitable for
279808	BUSTER Inflatable Collar, Nylon, XS	16 cm	5 cm	Yorkshire Terrier, Papillon etc.
279804	BUSTER Inflatable Collar, Nylon, S	22 cm	7 cm	Fox Terrier, Jack Russell, etc.
279805	BUSTER Inflatable Collar, Nylon, M	25.5 cm	11cm	Springer Spaniel, Border Collie etc.
279806	BUSTER Inflatable Collar, Nylon, L	30.5 cm	16 cm	Golden Retriever, Flat-coated etc.
279807	BUSTER Inflatable Collar, Nylon, XL	36 cm	18 cm	Rottweiler, Great Dane etc.
279809	BUSTER Inflatable Collar, Nylon, XXL	44 cm	22 cm	Old English Bulldog, Newfoundland etc.



KRUUSE Inflatable Collar PVC

Cat. No	Description	External diameter	Internal diameter	Suitable for
279814	KRUUSE Inflatable Collar, PVC, XS	14 cm	4 cm	Yorkshire Terrier, Papillon etc.
279815	KRUUSE Inflatable Collar, PVC, S	22 cm	7 cm	Fox Terrier, Jack Russell etc.
279816	KRUUSE Inflatable Collar, PVC, M	25.5 cm	11cm	Medium Spaniel, Border Collie etc.
279817	KRUUSE Inflatable Collar, PVC, L	30 cm	16 cm	Retriever, German Shepherd etc.
279818	KRUUSE Inflatable Collar, PVC, XL	44 cm	22 cm	Rottweiler, Great Dane etc.



BUSTER Body Suit & Body Sleeves for Dogs and Cats

The soft and comfortable BUSTER Body Suit and Body Sleeves are used as protection and comfort after surgery or in case of skin disease. The BUSTER Body Sleeves come in both a front leg and a hind leg version. It is an ideal solution for patients recovering from surgery, minor injuries to protect tender areas caused by excessive licking, biting or chewing.

BUSTER Body Suits are perfect for

- Postoperative care
- Wound care
- Spay/neuter recovery
- Hot spots and other skin ailments
- In cases of skin diseases
- Incontinence, urogenital or perianal issues
- In season

BUSTER Body Sleeves are perfect for

- Postoperative care
- Wound care
- Spay/neuter recovery
- Hot spots and other skin ailments





BUSTER Body Suit

The BUSTER Body Suit is available in a version for dogs and a version for cats.

Features and Benefits of the BUSTER Body Suit

- Comfortable, soft and flexible
- Freedom of movement
- Easy fitting and removal
- Press studs for more security
- Made from cotton/elastane

**DOUBLE FABRIC
UNDER THE ABDOMEN**



The BUSTER Body Suit is available in two versions for dogs and one version for cats.

BUSTER Body Suit EasyGo for Cats


- Double fabric under the abdomen enabling pads to be inserted in case of fluid loss
- Body Suit EasyGo for cats with extended tail opening
- The EasyGo version helps the cat when nature calls
- Suitable for female cats and tomcats

BUSTER Body Suit EasyGo for Dogs

- Body Suit EasyGo for dogs with extended tail opening and extra slit
- The EasyGo version helps the dog when nature calls
- Suitable for bitches and male dogs

BUSTER Body Suit Classic for Dogs

- Double fabric under the abdomen enabling pads to be inserted in case of fluid loss, incontinence or when the animal is in season
- Undo press studs and touch fasteners and pull back the fabric and roll up the fabric under the abdomen and above the back when nature calls
- Suitable for bitches and male dogs



BUSTER Body Suit


Instructions for Use

1. Insert the head into the BUSTER Body Suit
2. Insert each front leg separately into the openings
3. Pull the BUSTER Body Suit as far as possible over the body
4. Put the closing flap between the rear legs and insert the tail through the tail opening
5. Attach the closing flap on the touch fasteners on the loin and secure the press studs

Size Guide

The BUSTER Body Suit EasyGo for Cats is available in 4 sizes.


Size measured from collar to base of the tail.



Cat. No	Size	Body length	Chest circumference (2 cm from front legs)	Abdomen circumference (2 cm from hind legs)	Suitable for
273960	XXXS	27.5 cm	36 cm	30 cm	Young cat (small kitten)
273961	XXS	33 cm	42 cm	36 cm	Young cat (kitten)
273962	XS	38.5 cm	46 cm	40 cm	Adult cat
273963	S	43.5 cm	52 cm	46 cm	Adult cat (large breed)

The BUSTER Body Suit EasyGo for Dogs is available in 9 sizes.


Size measured from collar to base of the tail.



Cat. No	Size	Body length	Chest circumference (2 cm from front legs)	Abdomen circumference (2 cm from hind legs)	Suitable for
273951	XXXS	25 cm	30 cm	26 cm	Small Chihuahua
273952	XXS	32 cm	38 cm	34 cm	Chihuahua, Toy poodle
273953	XS	38 cm	46 cm	42 cm	Yorkshire, Terrier, Lhasa, Toy poodle
273954	S	41 cm	54 cm	50 cm	Malteser, Pekingeser, French Bulldog, Mops, Terrier
273955	M	45 cm	58 cm	52 cm	English Cocker Spaniel, Samojeed, Portuguese Podengo
273956	L	52 cm	62 cm	56 cm	Schnauzer, Beagle, Poodle, Fox Terrier, Chow-Chow
273957	XL	62 cm	70 cm	64 cm	Collie, Appenzeller, Boxer, Schnauzer, Labrador
273958	XXL	73 cm	78 cm	70 cm	Rottweiler, Golden Retriever, Herdershond
273959	XXXL	75 cm	86 cm	76 cm	Saint Bernhard, Newfoundlander

The BUSTER Body Suit Classic for Dogs is available in 9 sizes.

Size measured from collar to base of the tail.



Cat. No	Size	Body length	Chest circumference (2 cm from front legs)	Abdomen circumference (2 cm from hind legs)	Suitable for
273968	XXXS	26 cm	30 cm	26 cm	Small Chihuahua
273969	XXS	33 cm	38 cm	34 cm	Chihuahua, Toy poodle
273970	XS	39 cm	46 cm	42 cm	Yorkshire, Terrier, Lhasa, Toy poodle
273971	S	42 cm	54 cm	50 cm	Malteser, Pekingeser, French Bulldog, Mops, Terrier
273972	M	46 cm	58 cm	52 cm	English Cocker Spaniel, Samojeed, Portuguese Podengo
273973	L	53 cm	62 cm	56 cm	Schnauzer, Beagle, Poodle, Fox Terrier, Chow-Chow
273974	XL	63 cm	70 cm	64 cm	Collie, Appenzeller, Boxer, Schnauzer, Labrador
273975	XXL	74 cm	78 cm	70 cm	Rottweiler, Golden Retriever, Herdershond
273976	XXXL	76 cm	86 cm	76 cm	Saint Bernhard, Newfoundlander

The BUSTER Body Suit is washable at 40°C

BUSTER Body Sleeves Front Legs

BUSTER The Body Sleeves can be used for both cats and dogs. Please see the size guide.

Features and Benefits of the BUSTER Body Sleeves

- Comfortable, soft and flexible
- Breathable material
- Freedom of movement
- Easy fitting and removal
- Made from polyester/viscose/elastane
- Easy-to-cut non-fraying fabric




BUSTER Body Sleeves Front Legs

- The extra sleeve fabric can be folded around the paw of the dog in case of a paw injury. Please use an extra anti-slid paw protection for the safety of the dog.
- The remaining sleeve fabric can be cut off in a suitable length for the individual dog



BUSTER Body Sleeves Front Legs



Size Guide

The BUSTER Body Sleeves for front legs is available in 8 sizes.

Cat. No	Description	Sleeve Length	Body Length	Suitable for
274046	BUSTER Body Sleeves, front legs XXXS	21 cm	24 cm	Small Chihuahua, kitten
274032	BUSTER Body Sleeves, front legs, XXS	25 cm	30 cm	Chihuahua, Toy poodle, Shih Tzu, Cat
274033	BUSTER Body Sleeves, front legs, XS	29 cm	36 cm	Yorkshire Terrier, Lhasa Apso, Toy poodle, Cat
274034	BUSTER Body Sleeves, front legs, S	33 cm	42 cm	Maltese, Pekingese, French Bulldog, Pug, Terrier
274035	BUSTER Body Sleeves, front legs, M	37 cm	46 cm	English Cocker Spaniel, Samoyed, Portuguese Podengo
274036	BUSTER Body Sleeves, front legs, L	41 cm	50 cm	Schnauzer, Poodle, Fox Terrier, Chow-Chow
274037	BUSTER Body Sleeves, front legs, XL	45 cm	57 cm	Collie, Appenzeller, Boxer, Labrador Retriever
274038	BUSTER Body Sleeves, front legs, XXL	49 cm	64 cm	Rottweiler, Golden Retriever, Herdershond

Instructions for use

1. Roll up the BUSTER Body Sleeves between the hands having the KRUUSE logo at the top towards your own body
2. Pull the BUSTER Body Sleeves over the dog's head having the KRUUSE logo on top of the dog's body
3. Pull the BUSTER Body Sleeves over the dog's body
4. Roll up one of the sleeves with your hand from the bottom of the sleeve upwards
5. The dog's leg is lifted up and inserted into the sleeve. Your hand inside the sleeve holds the dog's paws avoiding the claws. to get caught. With your other hand you pull up the sleeve
6. Pull up the sleeve entirely till the elbow joint of the dog
7. Repeat bullet 4-6 on the second leg
8. Pull up the BUSTER Body Sleeves around the body
9. Put your hand along the dog's back inside the BUSTER Body Sleeves, from the back up to the neck, lift up enabling the BUSTER Body Sleeves to fit correctly into the elbow joint and under the body of the dog

It is important that the sleeves are pulled up in full so that the sleeves seam fully stretches on the dog's body and snuggle the dog. Adjust the length of the sleeves, the sleeves can be folded to suit the dog or the remaining fabric can be cut off.

The BUSTER Body Sleeves are machine-washable at 40°C

Instruction for use
Scan QR code to
watch the video



BUSTER Body Sleeves Hind Legs

The BUSTER Body Sleeves can be used for both cats and dogs. Please see the size guide.

Features and Benefits of the BUSTER Body Sleeves

- Comfortable, soft and flexible
- Breathable material
- Freedom of movement
- Easy fitting and removal
- Made from polyester/viscose/elastane
- Easy-to-cut non-fraying fabric




Let the dog keep the BUSTER Body Sleeves on when nature calls, and easily cut a customised slit or hole in the non-fraying fabric.

BUSTER Body Sleeves Hind Legs

- The extra sleeve fabric can be folded around the paw of the dog in case of a paw injury. Please use an extra anti-slid paw protection for the safety of the dog.
- The remaining sleeve fabric can be cut off in a suitable length for the individual dog



BUSTER Body Sleeves Hind Legs



Size Guide

The BUSTER Body Sleeves for hind legs is available in 8 sizes.

Cat. No	Description	Sleeve Length	Body Length	Suitable for
274047	BUSTER Body Sleeves, hind legs, XXXS	27 cm	28 cm	Small Chihuahua, kitten
274039	BUSTER Body Sleeves, hind legs, XXS	31 cm	36 cm	Chihuahua, Toy poodle, Shih Tzu, Cat
274040	BUSTER Body Sleeves, hind legs, XS	35 cm	45 cm	Yorkshire Terrier, Lhasa Apso, Toy poodle, Cat
274041	BUSTER Body Sleeves, hind legs, S	39 cm	50 cm	Maltese, Pekingese, French Bulldog, Pug, Terrier
274042	BUSTER Body Sleeves, hind legs, M	43 cm	55 cm	English Cocker Spaniel, Samoyed, Portuguese Podengo
274043	BUSTER Body Sleeves, hind legs, L	47 cm	61 cm	Schnauzer, Poodle, Fox Terrier, Chow-Chow
274044	BUSTER Body Sleeves, hind legs, XL	51 cm	71 cm	Collie, Appenzeller, Boxer, Labrador Retriever
274045	BUSTER Body Sleeves, hind legs, XXL	55 cm	81 cm	Rottweiler, Golden Retriever, Herdershond

Instructions for use

1. Roll up the BUSTER Body Sleeves and hold one sleeve between the hands, the KRUUSE logo must be on top of the dog's back
2. Put one of the dog's legs into the sleeve and pull the all the way up to till the dog's groin
Pull the tail through the tail hole at the end of the BUSTER Body Sleeves
3. Roll up the other sleeve and put the dog's second leg into the sleeve, pull up the sleeve entirely till the groin
4. Drag the BUSTER Body Sleeves over the body of the dog
5. Lift up one front leg and pull through the leg hole in the BUSTER Body Sleeves
6. Lift up the second front leg and pull it through the leg hole in the BUSTER Body Sleeves
7. Drag the BUSTER Body Sleeves all the way up to the neck of the dog. Put your hand inside the BUSTER Body Sleeves along the dog's back from the dog' ears to the tail - and pull the BUSTER Body Sleeves in place

It is important that the sleeves are pulled up in full so that the seam follows the body of the dog (in the groin).

Adjust the length of the sleeves of the hind legs, the sleeves can be folded to suit the dog or the remaining fabric can be cut off.

Instruction for use

Scan QR code to
watch the video



BUSTER Ear Cover

The soft and comfortable BUSTER Ear Cover provides both post-surgical and general head and ear protection. It can be used as protection in case of aural haematoma, cuts or tears and chronic itching or scratching.

The BUSTER Ear Cover is provided with strong, yet soft neoprene on the forehead with integrated flexible plastic spines for optimal fit. An extra hole can be cut in the easy-to-cut non-fraying fabric in order to leave out the healthy ear flap. The BUSTER Ear Cover is suitable for both cats and dogs.

Features and benefits

- Freedom of movement
- Easy fitting and removal
- Flexible and extremely breathable
- Extra straps for secure fastening to the dog's neck collar
- Soft and durable neoprene
- Machine-washable
- Available in 6 sizes and a starter pack

BUSTER Ear Cover is perfect for

- Recovery from aural hematoma
- Cuts or tears
- Chronic itching or scratching
- Bandage protection
- Wound care




Cut a hole in the non-fraying fabric for the healthy ear flap



Size Guideline

KRUUSE Cat. No	Description	Size	A Circumference around the forehead in cm	B Circumference around the neck in cm	C Length from forehead to neck in cm
161645	BUSTER Ear Cover, 1/pk	XS	20.5	24	12
161646	BUSTER Ear Cover, 1/pk	S	24	26	15
161647	BUSTER Ear Cover, 1/pk	M	27.5	29.5	17
161648	BUSTER Ear Cover, 1/pk	L	31.5	33.5	19.5
161649	BUSTER Ear Cover, 1/pk	XL	35	38	21.5
161650	BUSTER Ear Cover, 1/pk	2XL	38	42	24
161651	BUSTER Ear Cover, Starter Pack, 6/pk	XS - 2XL			



Material specification

Made from polyester/viscose/elastane and neoprene

The BUSTER Ear Cover is machine-washable at 40° C



Scan QR code to
watch the video



