






## KRUUSE Suture Needles - Regular Eye, 10/pk

Stainless steel suture needles in supreme quality regarding both strength and sharpness. Cutting needle is the preferred needle type when suturing in more fibrous tissues like the skin, whereas most surgeons prefer the round bodied needle when suturing soft tissues like intestines and bladder.





Regular eye

Cat. No	Shape	Length, mm	Symbol
152600	3/8 circle, reverse cutting	17 mm	 <p>3/8 Circle</p>
152601	3/8 circle, reverse cutting	21 mm	
152602	3/8 circle, reverse cutting	25 mm	
152603	3/8 circle, reverse cutting	29 mm	
152604	3/8 circle, reverse cutting	41 mm	
152605	3/8 circle, reverse cutting	49 mm	
152606	3/8 circle, reverse cutting	61 mm	
152607	3/8 circle, reverse cutting	89 mm (diam 1.1 mm)	
152608	3/8 circle, reverse cutting	89 mm (diam 1.6 mm)	
152620	1/2 circle, reverse cutting	41 mm	 <p>1/2 Circle</p>
152621	1/2 circle, reverse cutting	49 mm	
152622	1/2 circle, reverse cutting	69 mm	
152663	3/8 circle, conventional cutting	19 mm	 <p>3/8 Circle</p>
152664	3/8 circle, conventional cutting	25 mm	
152665	3/8 circle, conventional cutting	29 mm	
152666	3/8 circle, conventional cutting	39 mm	
152667	3/8 circle, conventional cutting	65 mm	
152668	3/8 circle, conventional cutting	89 mm	
152661	1/2 circle, conventional cutting	25 mm	
152662	1/2 circle, conventional cutting	39 mm	
152627	Straight, Cutting	60 mm	 <p>Straight</p>
152628	Straight, Cutting	90 mm	

KRUUSE Suture Needles - Regular Eye, 10/pk



Regular eye

Cat. No	Shape	Length, mm	Symbol
152625	3/8 Circle, Round Body, Taper point	28 mm	
152626	3/8 Circle, Round Body, Taper point	69 mm	
152610	1/2 Circle, Round Body, Taper point	23 mm	
152611	1/2 Circle, Round Body, Taper point	28 mm	
152612	1/2 Circle, Round Body, Taper point	38 mm	
152613	1/2 Circle, Round Body, Taper point	48 mm	
152614	1/2 Circle, Round Body, Taper point	53 mm	
152615	1/2 Circle, Round Body, Taper point	68 mm	
152617	1/2 Circle, Round Body, Taper cut	23 mm	