

# BUSTER Body Suit

The BUSTER Body Suit is available in a version for dogs and a version for cats.

## Features and Benefits of the BUSTER Body Suit

- Comfortable, soft and flexible
- Freedom of movement
- Easy fitting and removal
- Press studs for more security
- Made from cotton/elastane

**DOUBLE FABRIC  
UNDER THE ABDOMEN**



The BUSTER Body Suit is available in two versions for dogs and one version for cats.

### Body Suit EasyGo for Cats


- Double fabric under the abdomen enabling pads to be inserted in case of fluid loss.
- Body Suit EasyGo for cats with extended tail opening.
- The EasyGo version helps the cat when nature calls.
- It is suitable for female cats and tomcats.

### Body Suit EasyGo for Dogs

- Body Suit EasyGo for dogs with extended tail opening and extra slit.
- The EasyGo version helps the dog when nature calls.
- It is suitable for bitches and male dogs.

### Body Suit Classic for Dogs

- Double fabric under the abdomen enabling pads to be inserted in case of fluid loss, incontinence or when the animal is in season.
- Undo press studs and touch fasteners and pull back the fabric and roll up the fabric under the abdomen and above the back when nature calls.
- It is suitable for bitches and male dogs




# BUSTER Body Suit

## Instructions for Use

1. Insert the head into the Body Suit
2. Insert each front leg separately into the openings provided
3. Pull the Body Suit as far as possible over the body
4. Put the closing flap between the rear legs and insert the tail through the tail opening
5. Attach the closing flap on the touch fasteners on the loin and secure the press studs


## Guideline for Sizes

The BUSTER Body Suit EasyGo for cats is available in 4 sizes.  
Size measured from collar to base of the tail.




Cat. No	Size	Body length	Chest circumference (2 cm from front legs)	Abdomen circumference (2 cm from hind legs)	Suitable for
273960	3XS	27.5 cm	36 cm	30 cm	Young cat (small kitten)
273961	2XS	33 cm	42 cm	36 cm	Young cat (kitten)
273962	XS	38.5 cm	46 cm	40 cm	Adult cat
273963	S	43.5 cm	52 cm	46 cm	Adult cat (large breed)

The BUSTER Body Suit EasyGo for dogs is available in 9 sizes.  
Size measured from collar to base of the tail.



Cat. No	Size	Body length	Chest circumference (2 cm from front legs)	Abdomen circumference (2 cm from hind legs)	Suitable for
273951	3XS	25 cm	30 cm	26 cm	Small Chihuahua
273952	2XS	32 cm	38 cm	34 cm	Chihuahua, Toy poodle
273953	XS	38 cm	46 cm	42 cm	Yorkshire, Terrier, Lhasa, Toy poodle
273954	S	41 cm	54 cm	50 cm	Malteser, Pekingese, French Bulldog, Mops, Terrier
273955	M	45 cm	58 cm	52 cm	English Cocker Spaniel, Samojeed, Portuguese Podengo
273956	L	52 cm	62 cm	56 cm	Schnauzer, Beagle, Poodle, Fox Terrier, Chow-Chow
273957	XL	62 cm	70 cm	64 cm	Collie, Appenzeller, Boxer, Schnauzer, Labrador
273958	2XL	73 cm	78 cm	70 cm	Rottweiler, Golden Retriever, Herdershond
273959	3XL	75 cm	86 cm	76 cm	Saint Bernhard, Newfoundlander

The BUSTER Body Suit Classic for dogs is available in 9 sizes.  
Size measured from collar to base of the tail.



Cat. No	Size	Body length	Chest circumference (2 cm from front legs)	Abdomen circumference (2 cm from hind legs)	Suitable for
273968	3XS	26 cm	30 cm	26 cm	Small Chihuahua
273969	2XS	33 cm	38 cm	34 cm	Chihuahua, Toy poodle
273970	XS	39 cm	46 cm	42 cm	Yorkshire, Terrier, Lhasa, Toy poodle
273971	S	42 cm	54 cm	50 cm	Malteser, Pekingese, French Bulldog, Mops, Terrier
273972	M	46 cm	58 cm	52 cm	English Cocker Spaniel, Samojeed, Portuguese Podengo
273973	L	53 cm	62 cm	56 cm	Schnauzer, Beagle, Poodle, Fox Terrier, Chow-Chow
273974	XL	63 cm	70 cm	64 cm	Collie, Appenzeller, Boxer, Schnauzer, Labrador
273975	2XL	74 cm	78 cm	70 cm	Rottweiler, Golden Retriever, Herdershond
273976	3XL	76 cm	86 cm	76 cm	Saint Bernhard, Newfoundlander

The BUSTER Body Suit is washable at 40°C