

## BUSTER Ear Cover

The soft and comfortable BUSTER Ear Cover provides both post-surgical and general head and ear protection. It can be used as protection in case of aural haematoma, cuts or tears and chronic itching or scratching.

The BUSTER Ear Cover is provided with strong, yet soft neoprene on the forehead with integrated flexible plastic spines for optimal fit. An extra hole can be cut in the easy-to-cut non-fraying fabric in order to leave out the healthy ear flap. The BUSTER Ear Cover is suitable for both cats and dogs.

### Features and benefits

- Freedom of movement
- Easy fitting and removal
- Flexible and extremely breathable
- Extra straps for secure fastening to the dog's neck collar
- Soft and durable neoprene
- Machine-washable
- Available in 6 sizes and a starter pack

### BUSTER Ear Cover is perfect for

- Recovery from aural haematoma
- Cuts or tears
- Chronic itching or scratching
- Bandage protection
- Wound care




Cut a hole in the non-fraying fabric for the healthy ear flap



Scan QR code to watch the video

Size Guideline

KRUUSE Cat. No	Description	Size	A Circumference around the forehead in cm	B Circumference around the neck in cm	C Length from forehead to neck in cm
161645	BUSTER Ear Cover, 1/pk	XS	20.5	24	12
161646	BUSTER Ear Cover, 1/pk	S	24	26	15
161647	BUSTER Ear Cover, 1/pk	M	27.5	29.5	17
161648	BUSTER Ear Cover, 1/pk	L	31.5	33.5	19.5
161649	BUSTER Ear Cover, 1/pk	XL	35	38	21.5
161650	BUSTER Ear Cover, 1/pk	2XL	38	42	24
161651	BUSTER Ear Cover, Starter Pack, 6/pk	XS - 2XL			



Material specification

Made from polyester/viscose/elastane and neoprene

The BUSTER Ear Cover is machine-washable at 40° C

