

# BUSTER Premium Dog Collar

*The greatest innovation in e-collars since the BUSTER Collar!*



**Guideline for Sizes**

The BUSTER Premium Dog Collar is available in 5 sizes.

Cat. No	Packaging	Size	Body weight	Suitable for the following breeds
273528	5/pk	XX-small	Breeds from 1-3 kg	Chihuahua, Japanese Chin, Maltese, Pomeranian, Japanese Terrier, Yorkshire Terrier
273529	5/pk	X-small	Breeds from 3-6 kg	Bichon Frisé, Coton de Tulear, Papillon Dogs, American Hairless Terrier, Shih Tzu, Tibetan Spaniel, English Toy Spaniel
273531	5/pk	Small	Breeds 6-10 kg	Cairn Terrier, Miniature Schnauzer, Jack Russel Terrier, Dachshund, Poodle, Border Terrier, Cavalier King Charles Spaniel
273532	5/pk	Medium	Breeds 10-20 kg	Border Collie, Finnish Lapphund, Beagle, Staffordshire Bull Terrier, English Springer Spaniel, Cocker Spaniel
273533	5/pk	Large	Breeds 20-35 kg	German Shepherd, Labrador Retriever, Golden Retriever, Wirehaired Pointer, Boxer, Flatcoated Retriever
273534	5/pk		Starter Pack: 1 x XX-Small 1 x X-Small 1 x Small 1 x Medium 1 x Large	

**BUSTER Premium Dog Collar - Instructions for Use**



**Step 1**

- Insert the head into Dog Collar




**Step 2**

- Close the Velcro



*Easy handling - High comfort*



## BUSTER Premium Dog Collar

*The greatest innovation in e-collars since the BUSTER Collar!*

The soft and flexible BUSTER Premium Dog Collar is very comfortable and is used as protection after surgery or in case of skin disease.

The BUSTER Premium Dog Collar is an ideal choice for patients recovering from surgery or minor injuries. Also suitable to protect lesions resulting from or aggravated by excessive licking, biting or chewing.

The BUSTER Premium Dog Collar is perfect for

- Post operative care
- Wound care
- Spay/neuter recovery
- Hot spots and other skin lesions
- Skin disease

**WHEN ONLY THE BEST IS GOOD ENOUGH FOR THE PATIENT**

*Thought-through design ready to use*



**High comfort**  
The banding around the outer edge provides better comfort both for the patient and the surroundings



**Safety and perfect fit**  
The collar is provided with extra neoprene and Velcro closing for safety and perfect fit



**Easy to fasten**  
The collar has a safe and easy Velcro closing



The shield is made of 100% Polypropylene



**Avoid skin irritation**  
The banding around the inner edge provides better comfort and prevents skin irritation




**Comfortable perforated neoprene**  
The collar is provided with soft, flexible and comfortable perforated neoprene

## Dimensions of the BUSTER Premium Dog Collar

Size	XX-small	X-small	Small	Medium	Large
<b>A</b> Neck size					
Circumference	14 - 19 cm	19 - 25 cm	25 - 33 cm	33 - 38 cm	38 - 50 cm
Corresponding Ø (diameter)	4.5 - 6 cm	6 - 8 cm	8 - 10.5 cm	10.5 - 12 cm	12 - 16 cm
Neck size average <i>(All dimensions behind the ear)</i>	16.5 cm	22 cm	29 cm	35.5 cm	45.5 cm
<b>B</b> Shoulder size					
Circumference	22 cm	28 cm	37.5 cm	50 cm	63 cm
Corresponding Ø (diameter) <i>(Maximum expansion of collar at shoulder)</i>	7 cm	9 cm	12 cm	16 cm	20 cm
<b>C</b> Neck length <i>(Suitable span)</i>	3 - 7 cm	5 - 11 cm	8 - 14 cm	10 - 16.5 cm	14 - 22 cm
Average neck length	5 cm	8 cm	11 cm	14 cm	17 cm
<b>D</b> Length from ear to tip of the nose <i>(Suitable span)</i>	8 - 10 cm	10 - 13 cm	13 - 19 cm	19 - 25 cm	25 - 30 cm
Average length	9 cm	12 cm	16 cm	22 cm	27 cm

- A. Neck size
- B. Shoulder size
- C. Neck length
- D. Length from ear to tip of the nose



- **BETTER COMPLIANCE**
- **BETTER COMFORT**
- **READY TO USE**

