

BUSTER Easy Slide Cat Catheter



Popular cat catheter now available in BUSTER brand with the same good properties and then a little more...

Made from DEHP-free (phthalate free) medical grade of PTFE Polytetrafluoroethylene the BUSTER Easy Slide Cat catheter is available in two forms: with open end or side holes both suitable for drainage and/or bladder lavage.

Features and benefits

- High biocompatibility
- Suitable for medium-long placement (3 days)
- Optimal rigidity (no need for stylet), though still with proper softness
- Highly lubricious catheter tubing material that ensures easy catheterization
- Flexible baseplate with suturing holes allows securing of the catheter to the prepuce
- Luer Lock fitting which enables syringe-attachment for flushing and/or attaching a urine bag collection system to the catheter
- Radiopaque – visible in x-ray
- Individually packed in easy-open pouch
- Small package size

SMOOTH AND LUBRICIOUS FOR EASY AND UNCOMPLICATED CATHETERIZATION



Cat. No	Type	Size, Metric	Size, Imperial	Package
273400	Side holes	1.2 x 110 mm	3.5 Fr x 4 1/8"	5/pk
273401	Side holes	1.2 x 140 mm	3.5 Fr x 5 1/2"	5/pk
273436	Open end	1.2 x 110 mm	3.5 Fr x 4 1/8"	5/pk
273437	Open end	1.2 x 140 mm	3.5 Fr x 5 1/2"	5/pk
273438	Open end	1.2 x 180 mm	3.5 Fr x 7"	5/pk