


BUSTER Dog Catheter


The popular range of BUSTER Dog Catheter now consists of 7 sizes of catheters!

- Female luer lock fitting
- Rigid though flexible polypropylene catheter tubing for easy introduction
- Rounded end and atraumatic side holes to avoid tissue damage
- Userfriendly blister packaging
- 12/pk


Cat. No	Size, Metric	Size, Imperial
273405	1.3 mm x 50 cm	4 Fr x 20 in
273406	1.6 mm x 50 cm	5 Fr x 20 in
273410	2.0 mm x 50 cm	6 Fr x 20 in
273411	2.3 mm x 50 cm	7 Fr x 20 in
273420	2.6 mm x 50 cm	8 Fr x 20 in
273421	3.0 mm x 50 cm	9 Fr x 20 in
273430	3.3 mm x 50 cm	10Fr x 20 in



Rounded end and atraumatic side holes



Flexible polypropylene tubing



Female luer lock fitting

