

## BUSTER Barium Cat Catheter

The soft polyurethane catheter softens at body temperature and is less reactive compared to urine catheters made of other materials.

- Made of soft polyurethane and hence ideal for hospitalized patients
- Suitable for long term use (up to 3 days)
- Suture holes in the baseplate allows suturing to the prepuce
- Female luer lock fitting: easy connection to syringe or urine collection bag
- X-Ray detectable which makes the catheter ideal for locating obstructions
- Side holes and stylet
- 12/pk



Cat. No	Size, Metric	Size, Imperial
273467	1.0 x 130 mm	3 Fr x 5.2 in
273468	1.3 x 130 mm	4 Fr x 5.2 in