

BUSTER NC Urinary Catheter

-Non Coated

BUSTER NC Urinary Catheter is a soft and biocompatible catheter. It is made from DEHP- free (phthalate free) medical grade. The soft catheter is less likely to cause damages to the mucous membranes in urethra and in the bladder compared to the more rigid catheters.

- Atraumatic polished lateral drainage holes
- High biocompatibility
- Suitable for medium-long placement (3 days)
- Colour-coded luer lockfitting for easy size recognition
- User-friendly end-opening box
- 12/pk

FOR INCREASED PATIENT COMFORT

How to use...



Open the pouch



Add gel



The catheter is ready for catheterization

Assortment

Cat. No	Size, Metric	Size, Imperial	Colour
273875	1,6 mm x 50 cm	5 Fr x 20 in	Grey
273876	2,0 mm x 50 cm	6 Fr x 20 in	Green
273877	2,6 mm x 50 cm	8 Fr x 20 in	Blue
273878	3,3 mm x 50 cm	10 Fr x 20 in	Black

