

## BUSTER Urine bags 750 ml, Luer Lock Connection

Connect the urine bag to a BUSTER Catheter or BUSTER NC Urinary catheter to collect and control urine flow during catheterization.

The bag has a non-return valve, which prevent the urine from running retrograde. Further the bag has an indicative volume scale and a "push-pull" release valve, which makes it possible to empty the bag.

Luer Lock connection



Cat. No 273855  
10/pk

## KRUUSE Urine Bags, 750 ml, Adapter for Foley Catheter

Connect the urine bag to BUSTER Foley Catheter to collect and control urine flow during catheterization.

The bag has a non-return valve, which prevent the urine from running retrograde. Further the bag has an indicative volume scale and a "push-pull" release valve, which makes it possible to empty the bag.

Adapter for Foley Catheter



Cat. No 273856  
10/pk

## BUSTER Hanger for Urine Bags

Hanger for BUSTER Urine bags (Cat. No 273855 and 273856).

For easy handling of BUSTER Urine bags.

Easy to hang on surgery table or on the cage.

Cat. No 273857  
10/pk

