

# KRUUSE Postoperative Care





# KRUUSE Postoperative Care

This chapter contains our range of postoperative care products and rehabilitation equipment used by veterinarians and technicians worldwide. We offer a wide range of disposable collars for cats and dogs in the original BUSTER design and other products to prevent the animal from reaching certain areas of the body.

## BUSTER Collars

- BUSTER Premium Collar
- BUSTER Classic Collar
- BUSTER Classic Transparent Collar
- BUSTER Comfort Collar
- BUSTER Clic Collar
- BUSTER Design Collar



## BUSTER Premium Collar

The BUSTER Premium Dog Collar is the greatest innovation in e-collars since the original BUSTER Collar! Thought-through design ready to use.

The BUSTER Premium Dog Collar is an ideal choice for patients recovering from surgery or minor injuries. Also suitable to protect lesions resulting from or aggravated by excessive licking, biting or chewing.

### The BUSTER Premium Dog Collar is perfect for

- Post operative care
- Wound care
- Spay/neuter recovery
- Hot spots and other skin lesions
- Skin disease

### The BUSTER Premium Dog Collar is available in a Starter Pack containing

- 1 x size xx-small
- 1 x size x-small
- 1 x size small
- 1 x size medium
- 1 x size large

### Cat. No Description

273528	BUSTER Premium Dog Collar, xx-small, 5/pk
273529	BUSTER Premium Dog Collar, x-small, 5/pk
273531	BUSTER Premium Dog Collar, small, 5/pk
273532	BUSTER Premium Dog Collar, medium, 5/pk
273533	BUSTER Premium Dog Collar, large, 5/pk
273534	BUSTER Premium Dog Collar, Starter Pack



**BUSTER**

Cat. No	Packaging	Size	Body weight	Suitable for the following breeds
273528	5/pk	XX-Small	Breeds 1-3 kg	Chihuahua, Japanese Chin, Maltese
273529	5/pk	X-Small	Breeds 3-5 kg	Bichon Frisé, Coton de Tulear, Papillon Dogs
273531	5/pk	Small	Breeds 5-10 kg	Cairn Terrier, Poodle, Border Terrier, Cavalier King Charles Spaniel
273532	5/pk	Medium	Breeds 10 - 20 kg	Border Collie, Beagle, Staffordshire Bull Terrier, English Springer Spaniel
273533	5/pk	Large	Breeds 20 - 35 kg	German Shepherd, Labrador Retriever, Golden Retriever, Flatcoated Retriever

## BUSTER Clic Collar

Quick and simple clic closing system. Makes the application of the collar as easy as buttoning your shirt. Fully transparent material will allow the animal all-round vision. Soft and flexible.

Made of environmentally friendly material. 100% PP - polypropylene.

Cat. No	Description
273480	BUSTER Clic Collar 7.5 cm, 10/pk
273481	BUSTER Clic Collar 10 cm, 10/pk
273482	BUSTER Clic Collar 12.5 cm, 10/pk
273483	BUSTER Clic Collar 15 cm, 10/pk
273484	BUSTER Clic Collar 20 cm, 10/pk
273485	BUSTER Clic Collar 25 cm, 10/pk
273486	BUSTER Clic Collar 30 cm, 10/pk
273488	BUSTER Clic Collar 35 cm
273487	BUSTER Clic Collar 40 cm
273490	BUSTER combi pack of clic collar, 20/pk 2 x 7.5 cm, 3 x 10 cm, 3 x 12.5 cm, 3 x 15 cm, 3 x 20 cm, 3 x 25 cm and 3 x 30 cm
273489	BUSTER Clic Collar, assorted sizes from 7.5 to 30 cm, 7/pk



## BUSTER Comfort Collar

Collar with rubber edge - equipped with a new quick-fastening system. The collar is transparent which enables the animal to see through it. Soft and flexible.

Made of environmentally friendly material. 100% PP - polypropylene.

Cat. No	Description
273900	BUSTER Comfort Collar 7.5 cm, 5/pk
273901	BUSTER Comfort Collar 10 cm, 5/pk
273902	BUSTER Comfort Collar 12.5 cm, 5/pk
273903	BUSTER Comfort Collar 15 cm, 5/pk
273904	BUSTER Comfort Collar 20 cm, 5/pk
273905	BUSTER Comfort Collar 25 cm, 5/pk
273906	BUSTER Comfort Collar 30 cm, 5/pk
273907	BUSTER Comfort collars, assorted sizes from 7.5 to 30 cm, 7/pk




Illustration of BUSTER Clic / Comfort Collar



How to assemble the BUSTER Clic / Comfort Collar

# BUSTER Classic Collar, Transparent

Our well-known classic dog collar is also available in a transparent version. Flexible and easy to handle. Easy to fit on the patient; available in 9 sizes to ensure perfect fit for any dog.

Made of environmentally friendly material. 100% PP - polypropylene.

Cat. No	Description
273391	BUSTER Transparent Dog Collar 7.5 cm, 10/pk
273392	BUSTER Transparent Dog Collar 10 cm, 10/pk
273393	BUSTER Transparent Dog Collar 12.5 cm, 10/pk
273394	BUSTER Transparent Dog Collar 15 cm, 10/pk
273395	BUSTER Transparent Dog Collar 20 cm, 10/pk
273396	BUSTER Transparent Dog Collar 25 cm, 10/pk
273397	BUSTER Transparent Dog Collar 30 cm, 10/pk
273399	BUSTER Transparent Dog Collar 35 cm
273398	BUSTER Transparent Dog Collar 40 cm
273389	BUSTER Transparent Dog Collars, assorted sizes from 7.5 to 30 cm, 7/pk



# BUSTER Classic Collar

Protective disposable collars for dogs.

An inexpensive solution to prevent the patient from reaching and damaging wounds on head and body. Flexible and easy to handle. Easy to fit on the patient; available in 9 sizes to ensure perfect fit for any dog.

Made of environmentally friendly material. 100% PP - polypropylene.

Cat. No	Description
273310	BUSTER Classic Dog Collar 7.5 cm, 10/pk
273320	BUSTER Classic Dog Collar 10 cm, 10/pk
273325	BUSTER Classic Dog Collar 12.5 cm, 10/pk
273330	BUSTER Classic Dog Collar 15 cm, 10/pk
273340	BUSTER Classic Dog Collar 20 cm, 10/pk
273350	BUSTER Classic Dog Collar 25 cm, 10/pk
273360	BUSTER Classic Dog Collar 30 cm, 10/pk
273368	BUSTER Classic Dog Collar 35 cm
273367	BUSTER Classic Dog Collar 40 cm
273369	BUSTER Classic Dog Collars, assorted sizes from 7.5 to 30 cm, 7/pk





Illustration of BUSTER Classic Collar





How to assemble the BUSTER Classic Collar

# BUSTER Design Collar Bones

Make the sending home of the animal more fun with the new BUSTER Design Collar with printed bones. Watch the surprise and happiness of your clients when the dog wears printed collar in eye-catching colours.

- Quick BUSTER Clic Fastening system
- Color-coded according to size
- Full body protection
- Multiple adjustment points for proper fit
- 7 sizes available



## Recommended use for BUSTER Design Collar Bones

- Postoperative care
- Wound care
- Spay/neuter recovery
- Hot spots and other skin ailments
- Skin diseases




Cat. No	Size	Pk. Size	Colour
273731	7.5 cm	10/pk	
273732	10cm	10/pk	
273733	12.5 cm	10/pk	
273734	15 cm	10/pk	
273735	20 cm	10/pk	
273736	25 cm	10/pk	
273737	30 cm	10/pk	

## Instructions for use

Scan QR code to see the video



273538 BUSTER Design Collar Bones, assorted sizes from 7.5 to 30 cm, 7/pk



Size	Length in cm / Inch	Opening in cm / Inch	Neck-hole in cm / Inch	Suitable for	Kg
7.5 cm	7.5 cm / 2.95"	13 - 15 cm / 5.1" - 5.9"	7 - 8 cm / 2.7" - 3.1"	Pomeranian, Chihuahua	1-3 kg
10 cm	10 cm / 3.93"	15 - 18 cm / 5.9" - 7.0"	7 - 8 cm / 2.7" - 3.1"	Yorkshire terrier, Bichon Frise, Maltese	3-5 kg
12.5 cm	12.5 cm / 4.92"	19 - 21 cm / 7.4" - 8.2"	9 - 10.5 cm / 3.5" - 4.1"	Danish/Swedish farmdog, Papillon, Cairn Terrier, Border Terrier	5-10 kg
15 cm	15 cm / 5.90"	20 - 22 cm / 7.8" - 8.6"	10 - 12 cm / 3.9" - 4.7"	Dutch Kooikerhondje, Scottish Terrier, Beagle, Irish Terrier	10-15 kg
20 cm	20 cm / 8.87"	23 - 27 cm / 9.0" - 10.6"	12 - 14 cm / 4.7" - 5.5"	Australian Shepherd, English Springer Spaniel, Staffordshire Bull Terrier	15-25 kg
25 cm	25 cm / 9.84"	34 - 42 cm / 13.3" - 16.5"	14 - 16 cm / 5.5" - 6.2"	Labrador Retriever, Siberian Husky, Irish Setter	25-30 kg
30 cm	30 cm / 11.81"	41 - 48 cm / 16.1" - 18.8"	15 - 18 cm / 5.9" - 7.0"	Golden Retriever, Boxer, Collie, German Shepherd	25-35 kg

# BUSTER Design Collar Under Construction

The BUSTER Design Collar is an eye-catching printed collar with a roadworks theme and road signs with cute statements wishing the dog good luck.

The BUSTER Design Collar is transparent and has the quick and simple clic closing system, which is well known from the BUSTER Clic Collars.

Made of environmentally friendly material. 100% PP - polypropylene.

Cat. No	Description
273520	BUSTER Design Dog Collar 7.5 cm "Under Construction", 10/pk
273521	BUSTER Design Dog Collar 10 cm "Under Construction", 10/pk
273522	BUSTER Design Dog Collar 12.5 cm "Under Construction", 10/pk
273523	BUSTER Design Dog Collar 15 cm "Under Construction", 10/pk
273524	BUSTER Design Dog Collar 20 cm "Under Construction", 10/pk
273525	BUSTER Design Dog Collar 25 cm "Under Construction", 10/pk
273526	BUSTER Design Dog Collar 30 cm "Under Construction", 10/pk
273527	BUSTER Design Dog Collar "Under Construction", assorted sizes from 7.5 to 30 cm, 7/pk



# BUSTER Dog Collar for Brachycephalic Breeds

Our range of BUSTER Collars has been extended with the new BUSTER collar for brachycephalic breeds. Finally, you can get a suitable collar for the brachycephalic dogs without having to adjust and cut the collar.

- Quick fastening system
- Full body coverage
- Multiple adjustment points for proper fit

Recommended use for BUSTER dog collar for brachycephalic breeds:

- Postoperative care
- Wound care
- Spay/neuter recovery
- Hot spots and other skin ailments
- Skin diseases

The BUSTER dog collar for brachycephalic breeds is transparent and has the quick and simple "Clic" closing system which is well-known from the BUSTER Clic Collars.




The collars are made of 100% polypropylene.





Cat. No	Size	Pk. Size	Suitable for	Kg
273500	XS	10/pk	Pekingese, Japanese Chin, Brussels griffon, Chihuahua	2-6 kg
273501	S	10/pk	Pug, Boston terrier, Cavalier King Charles spaniel, Shih Tzu	6 – 9 kg
273502	M	10/pk	French Bulldog, Lhasa Apso	9 – 16 kg
273503	L	10/pk	English Bulldog, American bulldog, Boxer, Bullmastiff	18 - 25 kg

273504 BUSTER Collar for Brachycephalic Breeds, assorted sizes, 2 x XS, 3 x S, 3 x M, 2 x L



Size	Neck-hole cm	Neck-hole inches	Opening cm	Opening inches	Length cm	Length inches
XS	8 - 9 cm	3.1 - 3.5"	17 - 20 cm	6.6 - 7.8"	10 cm	3.9"
S	9.5 - 10.5 cm	3.7 - 4.1"	23 - 25 cm	9.0 - 9.8"	12 cm	4.7"
M	12 - 14 cm	4.7 - 5.5"	30 - 33 cm	11.8 - 12.9"	16 cm	6.2"
L	14 - 17 cm	5.5 - 6.6"	36 - 40 cm	14.1 - 15.7"	20 cm	7.8"

## Collar Rack

Stainless steel rack for storing and displaying the BUSTER collars in the clinic. Made for seven sizes.

Cat. No	Description
273370	Stainless steel rack for collars



## Strap to Fix Collars

To fix BUSTER collars or to be used as an identity band.

Cat. No	Description
273300	Strap to fix BUSTER collars, 10/pk



## BUSTER Neckstrap

Practical and flexible neckstrap, which is used in connection with a BUSTER collar when the dog's normally worn collar is not suitable for this purpose.

Cat. No	Description
273492	Neckstrap for 7.5 and 10 cm, BUSTER clic collar, 10/pk
273494	Neckstrap for 12.5 and 15 cm, BUSTER clic collar, 10/pk
273497	Neckstrap 20, 25 & 30 cm, BUSTER clic collar, 10/pk
273498	Neckstrap for 40 cm, BUSTER clic collar, 10/pk



## BUSTER Soft Collar

For cats and small dogs. Comfortable, soft, non-woven material which is waterrepellent and easy to clean. Particularly suitable for cats due to the shape and flexibility of the collar.

Cat. No	Description
273384	BUSTER Soft collar for cats and small dogs, 10/pk



## Rabbit Collar

Protective collar for rabbits. Made of flexible and strong polyethylene. An inexpensive solution to prevent the animals from reaching wounds on head and body.

Cat. No	Description
273390	BUSTER rabbit collar, 10/pk



## Bird Collar

Disposable collar for birds. Prevents the bird from picking feathers and wounds. Easy to fit on the bird with an office stapler. Transparent material.

Cat. No	Description
273375	BUSTER bird collars 80 mm, 10/pk
273376	BUSTER bird collars 100 mm, 10/pk
273377	BUSTER bird collars 140 mm, 10/pk
273378	BUSTER bird collars 200 mm, 10/pk
273372	BUSTER bird collars combi, 15/pk
273386	BUSTER parrot collar, normal, 35 mm
273387	BUSTER parrot collar, large, 45 mm

273386-87



273372-78



### Inflatable Collars

Soft and comfortable collar to prevent the dog from reaching injuries, stitches, rashes or wounds. The protective inflatable collar has an inflatable inner pouch that is inflated through an air valve. The inner pouch is protected by an outer nylon or PVC cover to resist bites and scratches.

#### Features

- Will not inhibit the animal's field of vision
- Durable long-term use
- Comfortable for the animal
- Protective collar against injuries, rashes and post-surgery procedures
- Will not mark or scrape furniture
- Scratch and bite-resistant



#### BUSTER Inflatable Collar Nylon

Available in 6 sizes:

Cat. No	Description	External diameter	Internal diameter	Suitable for
279808	BUSTER Inflatable Collar, Nylon, XS	16 cm	5 cm	Yorkshire Terrier, Papillon etc.
279804	BUSTER Inflatable Collar, Nylon, S	22 cm	7 cm	Fox Terrier, Jack Russell, etc.
279805	BUSTER Inflatable Collar, Nylon, M	25.5 cm	11cm	Springer Spaniel, Border Collie etc.
279806	BUSTER Inflatable Collar, Nylon, L	30.5 cm	16 cm	Golden Retriever, Flat-coated etc.
279807	BUSTER Inflatable Collar, Nylon, XL	36 cm	18 cm	Rottweiler, Great Dane etc.
279809	BUSTER Inflatable Collar, Nylon, XXL	44 cm	22 cm	Old English Bulldog, Newfoundland etc.



#### KRUUSE Inflatable Collar PVC

Available in 5 sizes:

Cat. No	Description	External diameter	Internal diameter	Suitable for
279814	KRUUSE Inflatable Collar, PVC, XS	14 cm	4 cm	Yorkshire Terrier, Papillon etc.
279815	KRUUSE Inflatable Collar, PVC, S	22 cm	7 cm	Fox Terrier, Jack Russell etc.
279816	KRUUSE Inflatable Collar, PVC, M	25.5 cm	11cm	Medium Spaniel, Border Collie etc.
279817	KRUUSE Inflatable Collar, PVC, L	30 cm	16 cm	Retriever, German Shepherd etc.
279818	KRUUSE Inflatable Collar, PVC, XL	44 cm	22 cm	Rottweiler, Great Dane etc.



# BUSTER Body Suit & Body Sleeves for Dogs & Cats

The soft and comfortable BUSTER Body Suit and Body Sleeves are used as protection and comfort after surgery or in case of skin disease. It is an ideal solution for patients recovering from surgery or minor injuries or have tender areas caused by excessive licking, biting or chewing.

## BUSTER Body Suits are perfect for

- Postoperative care
- Wound care
- Spay/neuter recovery
- Hot spots and other skin ailments
- In cases of skin diseases
- Incontinence, urogenital or perianal issues
- In season

## BUSTER Body Sleeves are perfect for

- Postoperative care
- Wound care
- Spay/neuter recovery
- Hot spots and other skin ailments





# BUSTER Body Suit

The BUSTER Body Suit is available in a version for dogs and a version for cats.

## Features and Benefits of the BUSTER Body Suit

- Comfortable, soft and flexible
- Freedom of movement
- Easy fitting and removal
- Press studs for more security
- Made from cotton/elastane



The BUSTER Body Suit is available in two versions for dogs and one version for cats.

### Body Suit EasyGo for Cats


- Double fabric under the abdomen enabling pads to be inserted in case of fluid loss
- Body Suit EasyGo for cats with extended tail opening
- The EasyGo version helps the cat when nature calls
- It is suitable for female cats and tomcats

### Body Suit EasyGo for Dogs

- Body Suit EasyGo for dogs with extended tail opening and extra slit
- The EasyGo version helps the dog when nature calls
- It is suitable for bitches and male dogs

### Body Suit Classic for Dogs

- Double fabric under the abdomen enabling pads to be inserted in case of fluid loss, incontinence or when the animal is in season
- Undo press studs and touch fasteners and pull back the fabric and roll up the fabric under the abdomen and above the back when nature calls
- It is suitable for bitches and male dogs





# BUSTER Body Suit


## Instructions for Use

1. Insert the head into the Body Suit
2. Insert each front leg separately into the openings provided
3. Pull the Body Suit as far as possible over the body
4. Put the closing flap between the rear legs and insert the tail through the tail opening
5. Attach the closing flap on the touch fasteners on the loin and secure the press studs

## Guideline for Sizes

The BUSTER Body Suit EasyGo for cats is available in 4 sizes.


Size measured from collar to base of the tail.



Cat. No	Size	Body length	Chest circumference (2 cm from front legs)	Abdomen circumference (2 cm from hind legs)	Suitable for
273960	3XS	27.5 cm	36 cm	30 cm	Young cat (small kitten)
273961	2XS	33 cm	42 cm	36 cm	Young cat (kitten)
273962	XS	38.5 cm	46 cm	40 cm	Adult cat
273963	S	43.5 cm	52 cm	46 cm	Adult cat (large breed)

The BUSTER Body Suit EasyGo for dogs is available in 9 sizes.


Size measured from collar to base of the tail.



Cat. No	Size	Body length	Chest circumference (2 cm from front legs)	Abdomen circumference (2 cm from hind legs)	Suitable for
273951	3XS	25 cm	30 cm	26 cm	Small Chihuahua
273952	2XS	32 cm	38 cm	34 cm	Chihuahua, Toy poodle
273953	XS	38 cm	46 cm	42 cm	Yorkshire, Terrier, Lhasa, Toy poodle
273954	S	41 cm	54 cm	50 cm	Malteser, Pekingeser, French Bulldog, Mops, Terrier
273955	M	45 cm	58 cm	52 cm	English Cocker Spaniel, Samojeed, Portuguese Podengo
273956	L	52 cm	62 cm	56 cm	Schnauzer, Beagle, Poodle, Fox Terrier, Chow-Chow
273957	XL	62 cm	70 cm	64 cm	Collie, Appenzeller, Boxer, Schnauzer, Labrador
273958	2XL	73 cm	78 cm	70 cm	Rottweiler, Golden Retriever, Herdershond
273959	3XL	75 cm	86 cm	76 cm	Saint Bernhard, Newfoundlander

The BUSTER Body Suit Classic for dogs is available in 9 sizes.

Size measured from collar to base of the tail.



Cat. No	Size	Body length	Chest circumference (2 cm from front legs)	Abdomen circumference (2 cm from hind legs)	Suitable for
273968	3XS	26 cm	30 cm	26 cm	Small Chihuahua
273969	2XS	33 cm	38 cm	34 cm	Chihuahua, Toy poodle
273970	XS	39 cm	46 cm	42 cm	Yorkshire, Terrier, Lhasa, Toy poodle
273971	S	42 cm	54 cm	50 cm	Malteser, Pekingeser, French Bulldog, Mops, Terrier
273972	M	46 cm	58 cm	52 cm	English Cocker Spaniel, Samojeed, Portuguese Podengo
273973	L	53 cm	62 cm	56 cm	Schnauzer, Beagle, Poodle, Fox Terrier, Chow-Chow
273974	XL	63 cm	70 cm	64 cm	Collie, Appenzeller, Boxer, Schnauzer, Labrador
273975	2XL	74 cm	78 cm	70 cm	Rottweiler, Golden Retriever, Herdershond
273976	3XL	76 cm	86 cm	76 cm	Saint Bernhard, Newfoundlander

The BUSTER Body Suit is washable at 40°C

## BUSTER Body Sleeves Front Legs

The unique Body Sleeves are available in a version for front legs and a version for hind legs depending on the area affected on the dog. The Body Sleeves can be used for both cats and dogs. Please see the guidelines for sizes.

### Features and Benefits of the BUSTER Body Sleeves

- Comfortable, soft and flexible
- Breathable material
- Freedom of movement
- Easy fitting and removal
- Made from polyester/viscose/elastane
- Easy-to-cut non-fraying fabric




### BUSTER Body Sleeves Front Legs

- The extra sleeve fabric can be folded around the paw of the dog in case of a paw injury. Please use an extra anti-slid paw protection for the safety of the dog.
- The remaining sleeve fabric can be cut off in a suitable length for the individual dog



# BUSTER Body Sleeves Front Legs



## Guideline for Sizes

The BUSTER Body Sleeves for front legs is available in 8 sizes.

Cat. No	Description	Sleeve Length	Body Length	Suitable for
274046	BUSTER Body Sleeves, front legs XXXS	21 cm	24 cm	Small Chihuahua, kitten
274032	BUSTER Body Sleeves, front legs, XXS	25 cm	30 cm	Chihuahua, Toy poodle, Shih Tzu, Cat
274033	BUSTER Body Sleeves, front legs, XS	29 cm	36 cm	Yorkshire Terrier, Lhasa Apso, Toy poodle, Cat
274034	BUSTER Body Sleeves, front legs, S	33 cm	42 cm	Maltese, Pekingese, French Bulldog, Pug, Terrier
274035	BUSTER Body Sleeves, front legs, M	37 cm	46 cm	English Cocker Spaniel, Samoyed, Portuguese Podengo
274036	BUSTER Body Sleeves, front legs, L	41 cm	50 cm	Schnauzer, Poodle, Fox Terrier, Chow-Chow
274037	BUSTER Body Sleeves, front legs, XL	45 cm	57 cm	Collie, Appenzeller, Boxer, Labrador Retriever
274038	BUSTER Body Sleeves, front legs, XXL	49 cm	64 cm	Rottweiler, Golden Retriever, Herdershond

## Instructions for use

1. Roll up the Body Sleeves between the hands having the KRUUSE logo at the top towards your own body
2. Pull the Body Sleeves over the dog's head having the KRUUSE logo on top of the dog's body
3. Pull the Body Sleeves over the dog's body
4. Roll up one of the sleeves with your hand from the bottom of the sleeve upwards
5. The dog's leg is lifted up and inserted into the sleeve. Your hand inside the sleeve holds the dog's paws avoiding the claws etc. to get caught. With your other hand you pull up the sleeve
6. Pull up the sleeve entirely till the elbow joint of the dog
7. Repeat bullet 4-6 on the second leg
8. Pull up the Body Sleeves around the body
9. Put your hand along the dog's back inside the Body Sleeve, from the back up to the neck, lift up enabling the Body Sleeve to fit correctly into the elbow joint and under the body of the dog

It is important that the sleeves are pulled up well so that the sleeve seam fully stretches the dog's body and snuggles the dog. Adjust the length of the sleeves, the sleeves can be folded to suit the dog or the remaining fabric can be cut off.

The BUSTER Body Sleeves are washable at 40°C



### Instructions for use

Scan QR code to see the video

# BUSTER Body Sleeves Hind Legs

The unique Body Sleeves are available in a version for front legs and a version for hind legs depending on the area affected on the dog.

The Body Sleeves can be used for both cats and dogs. Please see the guidelines for sizes.

## Features and Benefits of the BUSTER Body Sleeves

- Comfortable, soft and flexible
- Breathable material
- Freedom of movement
- Easy fitting and removal
- Made from polyester/viscose/elastane
- Easy-to-cut non-fraying fabric




Let the dog keep the Body Sleeves on when nature calls, and easily cut a customised slit or hole in the non-fraying fabric.

## BUSTER Body Sleeves Hind Legs

- The extra sleeve fabric can be folded around the paw of the dog in case of a paw injury. Please use an extra anti-slid paw protection for the safety of the dog.
- The remaining sleeve fabric can be cut off in a suitable length for the individual dog



# BUSTER Body Sleeves Hind Legs



## Guideline for Sizes

The BUSTER Body Sleeves for hind legs is available in 8 sizes.

Cat. No	Description	Sleeve Length	Body Length	Suitable for
274047	BUSTER Body Sleeves, hind legs, XXXS	27 cm	28 cm	Small Chihuahua, kitten
274039	BUSTER Body Sleeves, hind legs, XXS	31 cm	36 cm	Chihuahua, Toy poodle, Shih Tzu, Cat
274040	BUSTER Body Sleeves, hind legs, XS	35 cm	45 cm	Yorkshire Terrier, Lhasa Apso, Toy poodle, Cat
274041	BUSTER Body Sleeves, hind legs, S	39 cm	50 cm	Maltese, Pekingese, French Bulldog, Pug, Terrier
274042	BUSTER Body Sleeves, hind legs, M	43 cm	55 cm	English Cocker Spaniel, Samoyed, Portuguese Podengo
274043	BUSTER Body Sleeves, hind legs, L	47 cm	61 cm	Schnauzer, Poodle, Fox Terrier, Chow-Chow
274044	BUSTER Body Sleeves, hind legs, XL	51 cm	71 cm	Collie, Appenzeller, Boxer, Labrador Retriever
274045	BUSTER Body Sleeves, hind legs, XXL	55 cm	81 cm	Rottweiler, Golden Retriever, Herdershond

## Instructions for use

1. Roll up the Body Sleeve and one sleeve between the hands, the KRUUSE logo must be on top of the dog's back
2. Put one of the dog's leg legs into the sleeve and pull the sleeve entirely till the dog's groin  
Pull the tail through the tail hole at the end of the Body Sleeve
3. Roll up the other sleeve and put the dog's second leg into the sleeve, pull up the sleeve entirely till the groin
4. Drag the Body Sleeve over the body of the dog
5. Lift up one front leg and pull through the leg hole in the Body Sleeve
6. Lift up the second front leg and pull through the leg hole in the Body Sleeve
7. Drag the Body Sleeve all the way up to the neck of the dog. Put your hand inside the Body Sleeve along the dog's back from the dog's ears to the tail - and pull up the Body Sleeve in place

It is important that the sleeves are pulled up all the way so that the seam follows the body of the dog (in the groin).

Adjust the length of the sleeves of the hind legs, the sleeves can be folded to suit the dog or the remaining fabric can be cut off.



### Instructions for use

Scan QR code to see the video





