



LED NOW IN HEINE QUALITY.



HEINE G100 LED

Slit Illumination Head with LED Illumination

DURABILITY:

- > **HEINE offers a 5 year guarantee.**
HEINE makes no compromises in manufacturing high quality medical instruments. Our commitment to vertical integration in manufacturing means we control all aspects of our instrument quality, from utilising carefully selected, well-matched materials to a high level of hand assembly. This ensures that each HEINE instrument meets or exceeds all requirements in any medical environment.
- > **Solid metal construction and scratch-resistant swivel lens** guarantee a lifetime of maintenance-free use.
- > **LED with virtually unlimited hours of working life.** No need to ever replace bulbs.

PERFORMANCE:

- > **Slit Illumination Head.** Allows optimal instrumentation with magnification and illumination. Clears the way for access during examinations.
- > **Easy cleaning:** Round edges allow easy disinfection, also with standard disinfectants.
- > **Multi-coated, scratch-resistant swivel lens with 1.6x magnification** for a high-resolution image.
- > **Fiber Optic Illumination.** Ensures homogeneous, very bright illumination and an unobstructed view of the ear canal and tympanum.
- > **LED in HQ:** The new standard in LED illumination defining optimal light intensity, homogeneity and colour rendering for the most accurate diagnosis. Red is red, blue is blue. Colour temperature: Typ. 3,500K, Colour Rendering Index (CRI) >97, special index for red colours (R9) >93.
- > **Exclusive continuous brightness control** between 100% and 3% (patent pending) with one-finger operation.

FUNCTION:

- > **Double bayonet lock for specula.** Specula mounted firmly by double lock.
- > **For use in human and veterinary medicine.**
- > **The following specula and tips may be used:**
 - Long, closed metal specula
 - Slit metal specula, open
 - Spreadable speculum
 - Short, closed specula
 - UniSpec disposable specula

