

# KRUTEX Surgical Gowns

## Increased protection

- Ribbed cuffs
- Waterproof
- AAMI standards

## Comfortable fit

- Double side ties
- Durable fabric

## KRUTEX Premium Surgical Gown

- The KRUTEX Premium Surgical Gown is for use in procedures with a moderate risk of exposure
- Protection level AAMI 3\*
- The gown is made from non-woven SMMMS, a technology providing a high level of barrier protection and air permeability. It has a strong resistance to tears and punctures
- With the ultrasonic sealing the gown is free from threads meaning less material waste and no risk of thread deterioration
- Double internal and external side ties on the waist for optimal fit
- Ribbed cuff for increased protection
- Waterproof and surgically folded for ease of transportation and to preserve sterile surfaces
- Includes 2 hand towels (30 x 40 cm) for hand drying, cleanup, and for squaring off the surgical incision site, and to provide increased absorbency of blood and bodily fluids

Cat. No	Description
260446	KRUTEX Premium Surgical Gown, L=110 cm, S, 25/pk
260447	KRUTEX Premium Surgical Gown, L=125 cm, M, 25/pk
260448	KRUTEX Premium Surgical Gown, L=130 cm, L, 25/pk
260449	KRUTEX Premium Surgical Gown, L=145 cm, XL, 25/pk



## KRUTEX Basic Surgical Gown

- The KRUTEX Basic Surgical Gown is for use in procedures with a minimal to low risk of exposure
- Protection level AAMI 2\*
- The gown is made from non-woven SMS which features high tensile strength, softness and comfort for light duty applications
- Ribbed cuff for increased protection
- Waterproof and surgically folded for ease of transportation and to preserve sterile surfaces
- Includes 2 hand towels (10 x 40 cm) for hand drying, cleanup, and for squaring off the surgical incision site, and to provide increased absorbency of blood and bodily fluids

Cat. No	Description
260564	KRUTEX Basic Surgical Gown, L=110 cm, S, 50/pk
260565	KRUTEX Basic Surgical Gown, L=120 cm, M, 50/pk
260566	KRUTEX Basic Surgical Gown, L= 125 cm, L, 50/pk
260567	KRUTEX Basic Surgical Gown, L= 130 cm, XL, 50/pk



### \*AAMI Protection

AAMI provides standards for liquid barrier performance with ratings for different levels of gown protection from level 1-4.

**AAMI LEVEL 2** Minimal to low fluid barrier protection. For use for only low-fluid, minimally invasive surgical procedures, lumps and bumps.

**AAMI LEVEL 3** Moderate fluid barrier protection. Used for the widest range of surgical procedures, where moderate fluid protection is indicated.