

KRUUSE Autoclave Products

Optimise the daily autoclaving of instruments and trays

To cover most needs the KRUUSE autoclave product range consists of self-seal flat pouches, heat-seal rolls and gusseted pouch and autoclave tape. All products have changeable colour fields for easy verification of completed autoclave cycle.

KRUUSE Self-Seal Autoclave Pouches

The self-seal pouches make the daily handling of instruments and trays fast and easy, as there is no need for neither tape nor a welding machine to close the pouches.

In addition, all self-seal pouches have a writing area for easy application of sterilisation date.

| KRUUSE Cat. No | Size | Package |
|----------------|--------------|---------|
| 299250 | 90 x 165 mm | 200/pk |
| 299251 | 90 x 260 mm | 200/pk |
| 299252 | 100 x 400 mm | 200/pk |
| 299253 | 135 x 300 mm | 200/pk |
| 299254 | 180 x 340 mm | 200/pk |
| 299255 | 250 x 380 mm | 200/pk |
| 299256 | 250 x 500 mm | 200/pk |



KRUUSE Heat-Seal Gusseted Autoclave Pouch

If you have larger items or trays for autoclaving, we recommend that you choose the gusseted pouch as the height of the pouch makes insertion of the tray much easier.

Seal the gusseted pouch with a welding machine or tape.

| KRUUSE Cat. No | Size | Package |
|----------------|-------------------|---------|
| 299257 | 300 x 500 x 60 mm | 100/pk |



KRUUSE Autoclave Rolls

The range of rolls includes various sizes for both small and large instruments and trays. Seal the rolls with a welding machine or tape.

| KRUUSE Cat. No | Size | Package |
|----------------|----------------|---------|
| 299258 | 50 mm x 200 m | 1/pk |
| 299259 | 75 mm x 200 m | 1/pk |
| 299260 | 100 mm x 200 m | 1/pk |
| 299261 | 150 mm x 200 m | 1/pk |
| 299262 | 200 mm x 200 m | 1/pk |
| 299263 | 250 mm x 200 m | 1/pk |



KRUUSE Autoclave Tape

The autoclave tape has changeable colour fields for easy verification of completed autoclave cycle.

| KRUUSE Cat. No | Size | Package |
|----------------|--------------|---------|
| 299264 | 19 mm x 50 m | 1/pk |



Applicable to all KRUUSE Autoclave Products:

Temperature limit: 121°C for 20 minutes or 134°C for 5 minutes on steam.
Shelf life: 3 years

NB! After sterilisation the instruments will remain sterile for up to 6 months.