

KRUUSE Rehab Pro Knee Protector

All-in-one adaptable knee protector

Provides

- Compression
- Warmth
- Support

Suitable for

- Sore and painful knee joints
- Pre-surgery support for knee injuries
- Support for chronic knee problems
- Rehabilitation

Multiple use - one product covers it all

Useful for conditions like arthritis, arthrosis, patella luxation and traumatic knee injuries.

Soft Comfortable Support

Remove the side splints when soft support, warmth and compression are needed during rest time or for arthritis or conditions that include sore or stiff joints.

Rigid Support

Insert the side splints when rigid support is needed for acute, traumatic injuries, or when increased patella stability is necessary during walks or during restricted exercise programmes.

Perfect Fit

All straps are adjustable to offer a custom fit.

Consult the treating veterinarian for post-surgery use.



The overall design ensures extra comfort to the patient

The extra built-in features and soft fabric ensure that the dog is comfortable when wearing the KRUUSE Rehab Pro Knee Protector in daily activities, on walks, and when resting.



	XXS	XS	S	M	L	XL	XXL
A	28 cm	32 cm	36 cm	45 cm	53 cm	61.5 cm	71 cm
B	6 cm	8 cm	10 cm	12 cm	15 cm	17 cm	19 cm
C	9.5 cm	10.5 cm	11.5 cm	14.5 cm	16 cm	16.5 cm	19.5 cm
D	9 cm	10 cm	11 cm	14 cm	15 cm	16 cm	19.5 cm
E	4 cm	4 cm	4.5 cm	5.5 cm	6.5 cm	8 cm	9 cm

Size	Cat. No Left	Cat. No Right	Suitable for
XXS	279982	279974	Bichon Frisé, Coton de Tulear, Shih Tzu, Tibetan Spaniel
XS	279983	279975	Sheltie, Cocker Spaniel, Fox Terrier, Cavalier King Charles Spaniel
S	279984	279976	Kelpie, Whippet, Beagle, Boston Terrier
M	279985	279977	Labrador Retriever, Border Collie, Springer Spaniel, Australian Shepherd
L	279986	279978	German Shepherd, Female Rottweiler, Boxer, Siberian Husky
XL	279987	279979	Rottweiler, Bernese Mountain Dog, Akita, Bloodhound
XXL	279988	279980	Great Dane, Berner Sennen, Mastiff, Newfoundland, Irish Wolfhound