

KRUUSE range of Professional Stethoscopes

Our selection of professional stethoscopes are acoustically very sensitive to optimize auscultation

- All KRUUSE Stethoscopes have black tubing
- Each stethoscope includes a spare part kit with two sets of earplugs and one diaphragm



For Veterinary Use Only

KRUUSE Professional Stethoscope

- Designed for auscultation in general veterinary practise
- Acoustic sensitive
- Dual head
- Soft eartips
- Total length: 77 cm
- Large diaphragm: 44 mm
- Stainless steel chestpiece and binaural

Cat. No 292400



**KRUUSE Dual Tube Stethoscope
(also known as Rappaport or Sprague)**

- Designed for auscultation in general veterinary practise
- Acoustic sensitive
- Dual tubing
- Dual head
- Soft eartips
- Total length: 74.5 cm
- Large diaphragm: 40 mm
- Small diaphragm: 26 mm
- Zinc alloy chestpiece and aluminium binaural

Cat. No 292401



KRUUSE Small Head Stethoscope

- Designed to optimise auscultation in small animals
- Acoustic sensitive
- Dual small head
- Soft eartips
- Total length: 76.5 cm
- Large diaphragm: 40 mm
- Stainless steel chestpiece and binaural

Cat. No 292402



KRUUSE Cardiology Stethoscope

- Designed to improve auscultation in the low frequency area giving a better evaluation of heart sounds
- Acoustic sensitive
- Dual head
- Soft eartips
- Total length: 72.5 cm
- Large diaphragm: 44 mm
- Stainless steel chestpiece and aluminium binaural

Cat. No 292403

