

# Surgical Gowns

### Increased protection

- Ribbed cuffs
- Water-repellent
- AAMI standards

### Durable fabric

- High tensile strength
- Tear resistant
- Water-repellent

### Comfortable fit

- Double side ties
- Soft
- Breathable

## KRUUSE Reinforced Surgical Gown

Made from non-woven breathable laminated antistatic fabric (60 gsm) and SSMMS fabric (43 gsm) which ensures an optimal barrier.

- For procedures with a moderate risk of exposure
- Protection level AAMI 3\*
- Reinforced front for extra fluid protection
- Ultrasonic sealing avoids the risk of thread deterioration
- Ribbed cuffs for increased protection
- Water-repellent
- Surgically folded to preserve a sterile surface
- 2 hand towels (30 x 40 cm) for drying, absorbing, cleaning, and squaring off the surgical incision site



Cat. No	Description	Size	Package
281622	KRUUSE Reinforced Surgical Gown	M (L=125 cm)	25/pk
281623	KRUUSE Reinforced Surgical Gown	L (L=130 cm)	25/pk
281624	KRUUSE Reinforced Surgical Gown	XL (L=145 cm)	25/pk
281625	KRUUSE Reinforced Surgical Gown	XXL (L=159 cm)	25/pk



## KRUTEX Premium Surgical Gown

Made from non-woven SMMMS providing a high level of barrier protection and air permeability. Strong tear and puncture resistance.

- For procedures with a moderate risk of exposure
- Protection level AAMI 3\*
- Ultrasonic sealing avoids the risk of thread deterioration
- Double internal and external waist side ties for optimal fit
- Ribbed cuffs for increased protection
- Water-repellent
- Surgically folded to preserve a sterile surface
- 2 hand towels (30 x 40 cm) for drying, absorbing, cleaning, and squaring off the surgical incision site

Cat. No	Description	Size	Package
260446	KRUTEX Premium Surgical Gown	S (L=110 cm)	25/pk
260447	KRUTEX Premium Surgical Gown	M (L=125 cm)	25/pk
260448	KRUTEX Premium Surgical Gown	L (L=130 cm)	25/pk
260449	KRUTEX Premium Surgical Gown	XL (L=145 cm)	25/pk



## KRUTEX Basic Surgical Gown

Made from non-woven SMS which features high tensile strength, softness, and comfort for light duty applications.

- For procedures with a minimal to low risk of exposure
- Protection level AAMI 2\*
- Ribbed cuffs for increased protection
- Water-repellent
- Surgically folded to preserve a sterile surface
- 2 hand towels (10 x 40 cm) for drying, absorbing, cleaning, and squaring off the surgical incision site

Cat. No	Description	Size	Package
260564	KRUTEX Basic Surgical Gown	S (L=110 cm)	50/pk
260565	KRUTEX Basic Surgical Gown	M (L=120 cm)	50/pk
260566	KRUTEX Basic Surgical Gown	L (L=125 cm)	50/pk
260567	KRUTEX Basic Surgical Gown	XL (L=130 cm)	50/pk



### \*AAMI Protection

AAMI provides standards for liquid barrier performance with ratings for different levels of gown protection from level 1-4.



**AAMI LEVEL 2** Minimal to low fluid barrier protection. For use for only low-fluid, minimally invasive surgical procedures, lumps and bumps.



**AAMI LEVEL 3** Moderate fluid barrier protection. Used for the widest range of surgical procedures, where moderate fluid protection is indicated.