

KRUUSE Super-Absorb

Absorbent dressing with high absorbency properties for wounds with moderate to severe exudation. The excellent absorbent capacity reduces number of bandage changes and ensures better wound healing.

If the wound is highly infected, it is recommended to combine the dressing with systemic and/or topical antibiotic or antibacterial treatment.

Application

- Use as primary or secondary dressing for wounds of different sizes, shapes, and amount of exudate
- Effectively contains and encapsulates exudate and cellular debris to prevent micro-organisms from contaminating the wound
- The dressing is easily secured with adhesive tape or a bandage



Improved healing

- The 3 dressing layers contribute to improved healing without causing maceration
- The dressing protects the wound and absorbs wound exudate in a moist environment which promotes healing

More comfort

- Due to high absorption capacity and fluid retention, the dressing can stay on for days
- The dressing will not stick to the wound which minimises pain when it is removed or changed

Cat. No	Description	Size, Metric	Size, Imperial	Package
160319	KRUUSE Super-Absorb	10 x 10 cm	4" x 4"	10/pk
160320	KRUUSE Super-Absorb	10 x 20 cm	4" x 8"	10/pk
160321	KRUUSE Super-Absorb	20 x 20 cm	8" x 8"	10/pk
160322	KRUUSE Super-Absorb	20 x 30 cm	8" x 12"	10/pk



Scan QR code to watch the video

