

Absorbent Pads

KRUUSE absorbent pads effectively contains and encapsules fluids to ensure optimal patient comfort by keeping the patient dry. The pads are suitable for situations requiring hygienic and instant management of fluids.

Easy and safe usage

- Reduced cross contamination
- Easy to use and discard
- Improves working environment

Extra patient care

- Keeps the surface dry
- Soft and quilted material
- Very comfortable

KRUUSE Operation Pad

A unique non-sterile pad particularly suited for cleaning and as cage bedding. In 1 minute, the small pad absorbs 100 ml fluid, and the larger pad absorbs 200 ml.



Cat. No	Description	Size	Package
161282	KRUUSE Operating Pad	45 x 60 cm	20 x 10/pk
161283	KRUUSE Operating Pad	60 x 90 cm	12 x 10/pk

KRUUSE Dental Pad

A non-sterile absorbent pad particularly suited for dental examinations. The small pad absorbs up to 200 ml fluid in approx. 3 minutes and the larger pad absorbs up to 400 ml fluid in approx. 3-4 minutes.



Cat. No	Description	Size	Package
161248	KRUUSE Dental Pad	30 x 45 cm	20 x 10/pk
161249	KRUUSE Dental Pad	45 x 60 cm	20 x 10/pk