


KRUUSE Butterfly

- Designed to provide quick venous access

The IV Butterflies are ideal for blood sampling, injection or short term infusion

- Butterfly-shaped wings for easy handling
- Thin wall needles provide better flow rate per gauge
- Long bevel to minimize discomfort/pain during placement
- Tube with proximal, flexible luer lock fitting
- Colour-coded wings for instant identification of needle size
- Soft, non-toxic, non-irritant, medical grade tube
- Tube length 30 cm
- Sterile
- 50/pk

Cat. No	Item	Size	Size	Colour	Pk
121716	KRUUSE Butterfly	23 G x 3/4"	0.6 x 19 mm	Blue	50/pk
121717	KRUUSE Butterfly	21 G x 3/4"	0.8 x 19 mm	Green	50/pk
121718	KRUUSE Butterfly	19 G x 3/4"	1.1 x 19 mm	Yellow	50/pk



Long facet shaped bevel

