

KRUUSE InfuVein PRO IV Catheter

- Kink-resistant polyurethane tubing, that softens in body temperature
- Sharp needles with long bevel to minimize discomfort/pain during placement
- Transparent flash-back chamber
- Injection port (with coloured cap) provides easy opening and access to the vein
- The injection port's one way valve opens easily, and is ideal for administration of medication.
- We recommend use of 10 ml syringe or less with the injection port, otherwise the resistant will be too big
- The primary port can be used for fluid therapy
- Ideal for long term use (up to 7 days)
- Individually packed in blister
- 50/pk

Cat. No	Size	Gauge	Colour
121914	0.7 x 19 mm	24G	Yellow
121915	0.9 x 25 mm	22G	Blue
121916	1.1 x 32 mm	20G	Pink
121913	1.3 x 32 mm	18G	Green
121917	1.3 x 45 mm	18G	Teal

